

PHILIPS

Diagnostic Cardiology

Supplies catalog

Supplies catalog 2021





PHILIPS DigiTrak XT

PHILIPS
DigiTrak XT





Table of contents

Patient cables	2
Patient data cables	5
Replacement lead sets	6
Extended lead sets	7
Electrode connectors	7
Disposable ECG electrodes	8
Limb plate electrodes	9
Recording paper	10
Pouches and cases	11
Cardiograph batteries	12
List of host products	13

Patient cables

Ordering code	Description	Packaging	Use with	Europe USA
989803151631	12-Lead Complete Lead Set (AAMI), Standard Length (2 arm leads, 2 leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 2 Trident yoke labels (V1-V3, V4-V6), 4 colored clips, 14 colored rings, 6 brown base rings and IFU)	1 lead set	TC70: 860315, 860352; TC50: 860310, 860429; TC30: 860306	•
989803151641	12-Lead Complete Lead Set (IEC), Standard Length (2 arm leads, 2 leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 2 banana post adapters, 2 Trident yoke labels (C1-C3, C4-C6), 4 colored clips, 14 colored rings, 6 white base rings and IFU)	1 lead set		•
989803151651	12-Lead Complete Lead Set (AAMI), Long (2 long arm leads, 2 long leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 2 Trident yoke labels (V1-V3, V4-V6), 4 colored clips, 14 colored rings, 6 brown base rings and IFU)	1 lead set		•
989803151661	12-Lead Complete Lead Set (IEC), Long (2 long arm leads, 2 long leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 10 banana post adapters, 2 Trident yoke labels (C1-C3, C4-C6), 4 colored clips, 14 colored rings, 6 white base rings and IFU)	1 lead set		•
989803175891	12-Lead Patient Cable (IEC), Standard Length 3.4 m (11') long	1 patient cable	TC20: 860332	•
989803175911	12-Lead Patient Cable (IEC), Long 3.8 m (12.5') long	1 patient cable		•

Ordering code	Description	Packaging	Use with	Europe	USA
989803184921	12-Lead Patient Cable (IEC), Standard Length 3.4 m (11.2')	1 patient cable	TC10: 860392	•	
989803184941	12-Lead Patient Cable (IEC), Long 3.8 m (12.5') long	1 patient cable		•	
989803180121	ST80i 12-Lead Complete Lead Set (AAMI), Standard Length, Grabber	1 lead set	ST80i: 860343, 860351		•
989803180131	ST80i 12-Lead Complete Lead Set (IEC), Standard Length, Grabber	1 lead set		•	
989803180141	ST80i 12-Lead Complete Lead Set (AAMI), Standard Length, Snaps	1 lead set			•
989803180151	ST80i 12-Lead Complete Lead Set (IEC), Standard Length, Snaps	1 lead set		•	
989803180161	ST80i 12-Lead Complete Lead Set (AAMI), Long, Grabber	1 lead set			•
989803180171	ST80i 12-Lead Complete Lead Set (IEC), Long, Grabber	1 lead set		•	
989803180181	ST80i 12-Lead Complete Lead Set (AAMI), Long, Snaps	1 lead set			•
989803180191	ST80i 12-Lead Complete Lead Set (IEC), Long, Snaps	1 lead set		•	
989803157481	DigiTrak XT Patient Cable, Short	1 patient cable	DigiTrak XT: 860322	•	•
989803157491	DigiTrak XT Patient Cable, Medium	1 patient cable		•	•
989803157501	DigiTrak XT Patient Cable, Long	1 patient cable		•	•

Patient cables




989803151631.tif



989803175891.tif



989803175921



989803180151



989803184921.tif



989803184941.png



989803180191



989803157481
989803157491
989803157501

Patient data cables

Ordering code	Description	Packaging	Use with	Europe	USA
989803164281	USB Patient Data Cable Length: 2 m (6.6')	1 cable = 1 sales unit	TC70: 860315, 860316, 860352; TC50: 860310, 860311, 860429; TC30: 860306, 860307	•	•



989803164281.tif

Replacement lead sets

Ordering code	Description	Packaging	Use with	Europe	USA
989803151671	12-Lead Chest Lead Set (AAMI/IEC), Standard Length (1 x 1 into 3 Trident leads, 12 x 2 colored rings, 3 brown base rings, 3 white base rings (both AAMI and IEC), 3 banana post adapters, 4 Trident yoke lead description labels (V1-V3, V4-V6, C1-C3, C4-C6) and IFU)	1 lead set	TC70: 860315, 860352; TC50: 860310, 860429; TC30: 860306	•	•
989803151691	12-Lead Chest Lead Set (AAMI/IEC), Long (1 x 1 into 3 Trident leads, 3 x 2 colored rings, 3 brown base rings, 3 white base rings (both AAMI and IEC), 3 banana post adapters, 1 lead separator, 4 Trident yoke lead description labels (V1-V3, V4-V6, C1-C3, C4-C6) and IFU)	1 lead set		•	•
989803151711	12-Lead Limb Lead Set (AAMI/IEC), Standard Length (2 arm leads, 2 leg leads, 5 x 2 colored rings and 13 x 2 colored clips for AAMI and IEC, 4 banana post adapters and IFU)	1 lead set		•	•
989803151731	12-Lead Limb Lead Set (AAMI/IEC), Long (2 long arm leads, 2 long leg leads, 13 x 2 colored rings and 5 x 2 colored clips for AAMI and IEC, 4 banana post adapters and IFU)	1 lead set		•	•



989803151671.tif



989803151711.tif

Extended lead sets

Ordering code	Description	Packaging	Use with	Europe USA
989803151751	15/16/18-Lead Upgrade Kit (AAMI), Standard Length (4 additional lead assembly, 6 shorting plugs, 2 lead separators, 6 colored clips, 6 colored rings, 4 brown base rings and IFU)	1 set	TC70: 860315, 860352; TC50: 860310, 860429	•
989803151761	15/16/18-Lead Upgrade Kit (IEC), Standard Length (4 additional lead assembly, 6 shorting plugs, 2 lead separators, 4 banana post adapters, 6 colored clips, 6 colored rings, 4 white base rings and IFU)	1 set		•
989803151771	15/16/18-Lead Upgrade Kit (AAMI), Long (4 additional long lead assembly, 6 shorting plugs, 2 lead separators, 4 banana post adapters, 6 colored clips, 6 colored rings, 4 white base rings and IFU)	1 set		•
989803151781	15/16/18-Lead Upgrade Kit (IEC), Long (4 additional long lead assembly, 6 shorting plugs, 2 lead separators, 4 banana post adapters, 6 colored clips, 6 colored rings, 4 white base rings and IFU)	1 set		•

Electrode connectors

Ordering code	Description	Packaging	Use with	Europe USA
989803129241	Alligator clips, IEC For tab electrodes	10 adapters	TC70: 860315, 860352; TC50: 860310, 860429; TC30: 860306; TC20: 860332;	•
989803166031	Clear tab/snap adapter Enables the connection from 3–4 mm post banana plug coming from the cardiograph to tab and snap style electrodes	10 adapters	TC10: 860392	• •



989803166031.png

Disposable ECG electrodes

Ordering code	Alternative ordering code	Description	Packaging	Use with	Europe	USA
989803100441	13943B	Adult solid gel tab electrode for resting ECG 1-3/8 x 7/8" (35 x 22 mm)	100 per pack; 10 packs per case; = 1,000 electrodes	TC70: 860315, 860352; TC50: 860310, 860429; TC30: 860306; TC20: 860332; TC10: 860392	•	•
989803101301	40420A	High-tack foam, wet gel ECG electrode Used in general electrocardiographic procedures in particular patient ECG surveillance and ECG diagnosis recording. 1.1 x 2.2" (28 x 55 mm)	30 per pack; 10 packs per case; = 300 electrodes	ST80i: 860343, 860351 TC10: 860392	•	•
989803129091	M4612A	Solid gel snap ECG electrode	5 per pack; 60 packs per case; = 300 electrodes	DigiTrak XT: 860322	•	•
989803129101	M4613A	Solid gel snap electrode	30 per pack; 10 packs per case; = 300 electrodes		•	•
989803192551		Clear tape wet gel ECG electrode for monitoring adult ECG 50 mm (2") round	30 per pack; 10 packs per case; = 300 electrodes	DigiTrak XT: 860322	•	•



989803101301.jpg



989803129101



989803192551.jpg

Limb plate electrodes

Ordering code	Alternative ordering code	Description	Packaging	Use with	Europe USA
989803101691	40494E	Limb Clamp Electrode Clothespin type for 1/8" post end leads or 4 mm banana plug leads, screw connector	4 electrodes	All cardiographs with standard 1/8" post end or 4 mm banana plug leads	• •



989803101691.jpg

Recording paper

Ordering code	Alternative ordering code	Description	Packaging	Use with	Europe	USA
989803104771	M1708A	PageWriter XL thermal paper Z-fold with red grid and no header Size: Letter	1 case = 8 pads = 1,600 sheets	ST80i: 860343, 860351		•
989803104791	M1710A	PageWriter XL thermal paper Z-fold with red grid and no header Size: A4	1 case = 8 pads = 1,600 sheets		•	•
989803106261	M2481A	Thermal paper for PageWriter cardiographs Z-fold with 206 mm light red grid and 7.3 cm (2 7/8") header Size: Letter	1 case = 8 pads = 1,600 sheets	TC70: 860315, 860352; TC50: 860310, 860429; TC30: 860306; ST80i: 860343, 860351	•	•
989803106271	M2483A	Thermal paper for PageWriter cardiographs Z-fold with 233 mm light red grid and 6.7 cm (2.6") header Size: A4	1 case = 8 pads = 1,600 sheets		•	•
989803106281	M2485A	PTP brand anti-fade thermal paper for PageWriter cardiographs Z-fold with 206 mm light red grid and 7.3 cm (2 7/8") header Size: Letter	1 case = 8 pads = 1,600 sheets		•	•
989803136421	M3707A	Thermal paper for PageWriter Trim Z-fold with light red grid and 7.3 cm (2 7/8") header Size: Letter	1 case = 16 pads = 1,600 sheets		•	•
989803136431	M3708A	Thermal paper for PageWriter Trim Z-fold with light red grid and 7.3 cm (2 7/8") header Size: A4	1 case = 16 pads = 1,600 sheets		•	•
989803185271		Thermal roll paper, 110 mm (4.3") wide	1 case = 10 rolls	TC10: 860392	•	

Pouches and cases

Ordering code	Description	Packaging	Use with	Europe	USA
989803153451	DigiTrak XT pouch	1 pouch	DigiTrak XT: 860322	•	•
989803158191	DigiTrak XT belt clip	1 each		•	•
989803158201	DigiTrak XT belt clip – 10 pack	10 per pack		•	•
989803180801	Stress patient belt to secure patient interface module (PIM) for stress system (adjustable size)	1 each	ST80i: 860343, 860351	•	•
989803180811	Stress PIM holder	1 each	ST801: 860343, 860351, 453564277261, 453564277271	•	•
989803166501	PageWriter TC30 and TC50 Cardiograph keyboard membrane cover	5 per pack	TC50: 860310, 860429; TC30: 860306; TC20: 860332		
989803166511	PageWriter TC70 Cardiograph keyboard membrane cover	5 per pack	TC70: 860315, 860352		



989803153451.tif



989803166511.tif



989803180801.tif

Cardiograph batteries

Ordering code	Description	Packaging	Use with	Europe	USA
989803170371	Lithium-ion battery * Product available in China only	1 battery	TC70: 860315, 860352; TC50: 860310, 860429; TC30: 860306		
989803185291	PageWriter TC10 battery	1 battery	TC10: 860392	•	•
989803193431	Lithium ion battery 7.252 Wh	1 battery	ECG100	•	•
989803194541	Lithium-ion 9-cell battery pack	1 battery	TC70: 860315, 860352; TC50: 860310, 860429; TC30: 860306	•	•



989803185291.png



989803194541.png

Diagnostic Cardiology product list

Ordering code	Alternative ordering code	Description
860306		PageWriter TC30 cardiograph
860307		PageWriter TC30 upgrade
860310		PageWriter TC50 cardiograph
860311		PageWriter TC50 upgrade
860315		PageWriter TC70 cardiograph
860316		PageWriter TC70 upgrade
860322		DigiTrak XT Holter recorder
860332		PageWriter TC20 cardiograph
860333		PageWriter TC20 upgrade
860343		ST80i Stress Test System
860351		ST80i Stress Test System upgrade
860352		PageWriter TC70 government bundle with trolley
860392		PageWriter TC10 cardiograph
860393		PageWriter TC10 upgrade
860429		PageWriter TC50 government bundle with trolley

