

PHILIPS

Patient Monitoring

Supplies catalog 2023



Table of contents

| | | | |
|-----------|------------------------------|------------|---|
| 2 | Gas monitoring | 61 | Non-invasive blood pressure |
| 2 | Anesthesia gas | 61 | Multi-patient cuffs |
| 4 | Spirometry CO ₂ | 68 | Single-patient cuffs |
| 5 | Transcutaneous gas | 75 | Accessories and hoses |
| 6 | Capnography | 76 | Invasive blood pressure |
| 6 | Sidestream EtCO ₂ | 78 | Paper |
| 18 | Mainstream EtCO ₂ | 79 | Pulse oximetry |
| 24 | Cardiac output | 79 | Reusable sensors |
| 26 | ECG | 82 | Single-patient sensors |
| 26 | Electrodes | 84 | Adapter cables |
| 36 | Cables and leads | 90 | IntelliVue Cableless measurement |
| 56 | Telemetry adapters | 92 | Temperature |
| 57 | EEG | 92 | Continuous measurement |
| 60 | NMT | 95 | Spot check measurement |
| | | 96 | Batteries |
| | | 98 | Accessories |
| | | 100 | Patient monitors and defibrillators product list |

Anesthesia gas monitoring

Anesthesia gas

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|----------------|--|--------|-----|
| 451261001391 | M1662A | Anesthetic cal gas (use with M1659A) Volume of gas: 6 liters Shelf life: 2 years from the date of manufacture | 1 cylinder | Anesthesia gas module: 862286, 862119; G1: 862165; G5: 862227 | • | • |
| 989803100361 | 13907A | Calibration tube assembly | 1 assembly | Anesthesia gas analyzer: 453563372301, 453563372311; Microstream [®] ext.: 862393, 865377 | • | • |
| 989803100841 | 15210-64010 | Gas: Cal 1 cylinders for tcpCO ₂ | 6 cylinders | Microstream [®] ext.: 862393, 865377 | • | • |
| 989803100851 | 15210-64020 | Gas: Cal 2 cylinders for tcpCO ₂ | 6 cylinders | | • | • |
| 989803104681 | M1659A | Calibration bag (use with M1660A and M1662A) | 1 bag | Anesthesia gas module: 862119; G1: 862165; G5: 862227 | • | • |
| 989803110871 | M1657B | Water trap (use with M1658A) | 25 water traps | G1: 862165; G5: 862227 | • | • |
| 989803191081 | | Water trap | 12 water traps | G7m: 866173 | • | • |

Anesthesia gas monitoring

Anesthesia gas



451261001391
M1662A

989803100361
13907A

989803100841
15210-64010

989803100851
15210-64020

989803104681
M1659A

989803110871
M1657B

989803191081

Spirometry CO₂ IntelliVue Spirometry CO₂ sensors

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|------------|---------------------------|--------|-----|
| 989803144661 | M2785A | Flow sensor, adult/pediatric, single-patient use ET: > 4.0 mm (0.16 in) | 10 sensors | Spirometry module: 865021 | • | • |
| 989803144671 | M2786A | Flow sensor, infant/neonatal, single-patient use ET: > 4.0 mm (0.16 in) | 10 sensors | | • | • |



989803144661
M2785A



989803144671
M2786A

Transcutaneous gas monitoring Transcutaneous gas

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-------------|---|--------|-----|
| 989803100841 | 15210-64010 | Cal 1 gas cylinders for calibration (meets international shipping requirements) | 6 cylinders | tcpO ₂ module: 862283; Calibration unit: 15210B; tcpO ₂ transducer: Microstream® ext.: 862393, 865377 | • | • |
| 989803100851 | 15210-64020 | Cal 2 gas cylinders for calibration (meets international shipping requirements) | 6 cylinders | | • | • |



989803100841
15210-64010



989803100851
15210-64020

Capnography

Sidestream CO₂ LoFlo Sidestream CO₂ sensor

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|-----------------------------------|-----------|--|--------|-----|
| 989803144591 | M2741A | Sidestream CO ₂ sensor | 1 sensor | Capnography ext.: 862187; MP2: 865040; MP5: 865024; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; HeartStart XL+: 861290 | • | • |



989803144591
M2741A

Capnography

Sidestream CO₂ LoFlo Sidestream CO₂, intubated

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|-------------------------------------|-----------|---|--------|-----|
| 989803206711 | n/a | Airway Adapter - LoFlo - ET >4.0 mm | 10 sets | Sidestream LoFlo sensor: M2741A with Capnography ext.: 862187; MP2: 865040; MP5: 865024; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; HeartStart XL+: 861290 | • | • |
| 989803206721 | n/a | Airway Adapter - LoFlo - ET ≤4.0 mm | 10 sets | | • | • |



989803206711



989803206721

Capnography

Sidestream CO₂
LoFlo Sidestream CO₂, intubated with humidification

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|-----------------------------------|-----------|---|--------|-----|
| 989803206731 | n/a | Airway Adapter Set H - ET >4.0 mm | 10 sets | Sidestream LoFlo sensor: M2741A with Capnography ext.: 862187; MP2: 865040; MP5: 865024; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; HeartStart XL+: 861290 | • | • |
| 989803206741 | n/a | Airway Adapter Set H - ET ≤4.0 mm | 10 sets | | • | • |



989803206731



989803206741

Capnography

Sidestream CO₂
LoFlo Sidestream CO₂, non-intubated with O₂ delivery

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-------------|---|--------|-----|
| 989803206641 | n/a | Cannula - LoFlo - CO ₂ /O ₂ Nasal - Adult | 10 cannulas | Sidestream LoFlo sensor: M2741A with Capnography ext.: 862187, X2: 865039; MP2: 865040; MP5: 865024 | • | • |
| 989803206651 | n/a | Cannula - LoFlo - CO ₂ /O ₂ Nasal - Pediatric | 10 cannulas | | • | • |
| 989803206661 | n/a | Cannula - LoFlo - CO ₂ /O ₂ Nasal - Infant | 10 cannulas | | • | • |
| 989803206691 | n/a | Cannula - LoFlo - CO ₂ /O ₂ Oral/Nasal-Adult | 10 cannulas | | • | • |
| 989803206701 | n/a | Cannula - LoFlo - CO ₂ /O ₂ Oral/Nasal-Pediatric | 10 cannulas | | • | • |



989803206641



989803206651



989803206661



989803206691

Capnography

Sidestream CO₂ LoFlo Sidestream CO₂, non-intubated without O₂ delivery

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-------------|--|--------|-----|
| 989803206611 | n/a | Cannula - LoFlo - CO ₂ Nasal - Adult | 10 cannulas | Sidestream LoFlo sensor:: M2741A with Capnography ext.: 862187, X2: 865039; MP2: 865040; MP5: 865024 | • | • |
| 989803206621 | n/a | Cannula - LoFlo - CO ₂ Nasal - Pediatric | 10 cannulas | | • | • |
| 989803206631 | n/a | Cannula - LoFlo - CO ₂ Nasal - Infant | 10 cannulas | | • | • |
| 989803206671 | n/a | Cannula - LoFlo - CO ₂ Oral/Nasal - Adult | 10 cannulas | | • | • |
| 989803206681 | n/a | Cannula - LoFlo - CO ₂ Oral/Nasal-Pediatric | 10 cannulas | | • | • |



989803206611



989803206621



989803206631



989803206671



989803206681



Capnography

Sidestream CO₂
Sidestream CO₂ supplies

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|--|--------|-----|
| 989803206751 | n/a | Sample Line - LoFlo - w/ Male Luer | 10 lines | Sidestream LoFlo sensor: M2741A with Capnography ext.: 862187; MP2: 865040; MP5: 865024; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; HeartStart XL+: 861290 | • | • |
| 989803206761 | n/a | Sample Line- LoFlo- w/Male Luer-w/Nafion | | | • | • |
| 989803199061 | n/a | LoFlo extension line, length: 1.52 m (60 in) | 25 lines | Sidestream LoFlo sensor:: M2741A with Capnography ext.: 862187; MP2: 865040; MP5: 865024; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; HeartStart XL+: 861290 | • | |



989803206751



989803206761



989803199061

Capnography

Sidestream CO₂
CapnoTrak sidestream CO₂ supplies

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|--|--------|-----|
| 989803199011 | n/a | Water filter assembly | 10 | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; Heartstart Intrepid: 867172 | • | |
| 989803199021 | n/a | CO ₂ sampling extension line | 25 lines | | • | |
| 989803199031 | n/a | O ₂ delivery extension line | 25 lines | | • | |
| 989803199041 | n/a | Dehumidification tubing | 10 tubes | | • | |
| 989803199061 | n/a | LoFlo extension line, length: 1.52 m (60 in) | 25 lines | | • | |



989803199011



989803199021



989803199031



989803199041

Capnography

Sidestream CO₂
CapnoTrak Sidestream CO₂, intubated, short-term

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|--|--------|-----|
| 989803198991 | n/a | CapnoTrak Sidestream CO ₂ airway adapter set, intubated, short-term ET: > 4.0 mm (0.16 in) | 25 sets | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; Heartstart Intrepid: 867172 | • | |
| 989803199001 | n/a | CapnoTrak Sidestream CO ₂ airway adapter set, intubated, short-term ET: <= 4.0 mm (0.16 in) | 25 sets | | • | |



989803198991



989803199001

Capnography

Sidestream CO₂
CapnoTrak Sidestream CO₂, non-intubated, short-term, with O₂ delivery

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-------------|--|--------|-----|
| 989803198921 | n/a | CO ₂ / O ₂ nasal cannula, large | 25 cannulas | CM10: 863301, 863365; CM12: 863303, 863366; | • | |
| 989803198931 | n/a | CO ₂ / O ₂ nasal cannula, medium | 25 cannulas | CM100: 863300; CM120: 863302; CM150: 863304; | • | |
| 989803198941 | n/a | CO ₂ / O ₂ nasal cannula, small | 25 cannulas | Heartstart Intrepid: 867172 | • | |
| 989803198971 | n/a | CO ₂ / O ₂ oral-nasal cannula, large | 25 cannulas | | • | |
| 989803198981 | n/a | CO ₂ / O ₂ oral-nasal cannula, medium | 25 cannulas | | • | |



989803198921



989803198931



989803198941



989803198971



989803198981

Capnography

Sidestream CO₂

CapnoTrak Sidestream CO₂, non-intubated, short-term, without O₂ delivery

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-------------|--|--------|-----|
| 989803198891 | n/a | CO ₂ nasal cannula, large | 25 cannulas | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; Heartstart Intrepid: 867172 | • | |
| 989803198901 | n/a | CO ₂ nasal cannula, medium | 25 cannulas | | • | |
| 989803198911 | n/a | CO ₂ nasal cannula, small | 25 cannulas | | • | |
| 989803198951 | n/a | CO ₂ oral-nasal cannula, large | 25 cannulas | | • | |
| 989803198961 | n/a | CO ₂ oral-nasal cannula, medium | 25 cannulas | | • | |



989803198891



989803198901



989803198911



989803198951



989803198961



Short term use: procedure/emergency

Microstream™ Advance sampling lines

| Ordering code | Description | Packaging | CE | FDA |
|------------------------------|--|--|------------|-----|
| Non-intubated patient | | | | |
| 989803204571 | Oral/nasal sampling line; samples from both the mouth and the nose; integrated O ₂ tubing for supplemental oxygen delivery | Adult; O ₂ male connector; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204371 | | Adult; O ₂ male connector; long: 4 m (13 ft) | 25 per box | • • |
| 989803204551 | | Adult; O ₂ tubing female connector; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204361 | | Adult; O ₂ tubing female connector; long: 4 m (13 ft) | 25 per box | • • |
| 989803204561 | Oral/nasal sampling line; samples from both the mouth and the nose; available with or without O ₂ delivery | Pediatric; no O ₂ delivery; standard: 2 m (6.5 ft) no O ₂ delivery; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204351 | | Pediatric; O ₂ tubing female connector; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204351 | | Pediatric; O ₂ tubing female connector; long: 4 m (13 ft) | 25 per box | • • |
| 989803204481 | Nasal sampling line; nasal sampling only; integrated O ₂ tubing for supplemental oxygen delivery | Adult; O ₂ tubing female connector; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204491 | | Adult; O ₂ tubing female connector; long: 4 m (13 ft) | 25 per box | • • |
| 989803204501 | | Pediatric; O ₂ tubing female connector; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204611 | Nasal sampling line; nasal sampling only; no O ₂ delivery | Adult; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204441 | | Infant/neonate; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204601 | | Pediatric; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204451 | 60fr bite block sampling line; samples from both the mouth and the nose; integrated O ₂ tubing for supplemental oxygen delivery; O ₂ delivery system can administer oxygen flow up to 10 L/min | O ₂ male connector; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204461 | | O ₂ tubing female connector; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204471 | | O ₂ tubing female connector; long: 4 m (13 ft) | 25 per box | • • |
| Intubated patient | | | | |
| 989803204511 | Intubated sampling line; for use with ETT size > 4.5 mm | Adult/pediatric; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204321 | | Adult/pediatric; long: 4 m (13 ft) | 25 per box | • • |

Microstream™ Advance sampling lines




Extended duration

Microstream™ Advance sampling lines

| Ordering code | Description | Packaging | CE | FDA |
|------------------------------|--|---|------------|-----|
| Non-intubated patient | | | | |
| 989803204381 | Oral/nasal sampling line; samples from both the mouth and the nose; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery | Adult; O ₂ tubing female connector; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204391 | Oral/nasal sampling line; samples from both the mouth and the nose; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery | Adult; O ₂ tubing female connector; long: 4 m (13 ft) | 25 per box | • • |
| 989803204401 | Oral/nasal sampling line; samples from both the mouth and the nose; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery | Pediatric; O ₂ tubing female connector; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204411 | Oral/nasal sampling line; samples from both the mouth and the nose; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery | Pediatric; O ₂ tubing female connector; long: 4 m (13 ft) | 25 per box | • • |
| 989803204621 | Nasal sampling line; nasal sampling only; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery available | Adult; no O ₂ delivery; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204581 | Nasal sampling line; nasal sampling only; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery available | Adult; O ₂ tubing female connector; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204591 | Nasal sampling line; nasal sampling only; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery available | Pediatric; O ₂ tubing female connector; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204631 | Nasal sampling line; nasal sampling only; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery available | Neonatal/infant; no O ₂ delivery; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204421 | Nasal sampling line; nasal sampling only; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery available | Neonatal/infant; O ₂ tubing female connector; medium: 3 m (9.8 ft) | 25 per box | • • |
| 989803204431 | Nasal sampling line; nasal sampling only; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery available | Neonatal/infant; no O ₂ delivery; long: 4 m (13 ft) | 25 per box | • • |
| Intubated patient | | | | |
| 989803204521 | Intubated sampling line; for use with ETT > 4.5 mm; drier moisture reduction technology | Adult/pediatric; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204331 | Intubated sampling line; for use with ETT > 4.5 mm; drier moisture reduction technology | Adult/pediatric; long: 4 m (13 ft) | 25 per box | • • |
| 989803204531 | Intubated sampling line; for use with ETT ≤ 4.5 mm; drier moisture reduction technology | Neonatal/infant; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204341 | Intubated sampling line; for use with ETT ≤ 4.5 mm; drier moisture reduction technology | Neonatal/infant; long: 4 m (13 ft) | 25 per box | • • |
| 989803204301 | Intubated sampling line; for use with ETT > 4.5 mm; longer drier moisture reduction technology for high ambient humidity | Adult/pediatric; standard: 2 m (6.5 ft) | 25 per box | • • |
| 989803204311 | Intubated sampling line; for use with ETT ≤ 4.5 mm; longer drier moisture reduction technology for high ambient humidity | Neonatal/infant; standard: 2 m (6.5 ft) | 25 per box | • • |

Microstream™ Advance sampling lines



Capnography

Mainstream CO₂ solution

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|----------------|--|--------|-----|
| 989803142651 | n/a | Mainstream CO ₂ sensor | 1 sensor | Capnography ext.: 862187; MP2: 865040; MP5: 865024, 865322; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; HeartStart XL+: 861290 | • | • |
| 989803142661 | n/a | Single-patient use adult airway adapter ET: > 4.0 mm (0.16 in) | 10 adapters | | • | • |
| 989803142671 | n/a | Single-patient use infant airway adapter ET: < 4.0 mm (0.16 in) Deadspace: < 1 cc | | | • | • |
| 989803206591 | n/a | Airway Adapter-Reusable-Adult/Pediatric | 1 adapter | | • | • |
| 989803206601 | n/a | Airway Adapter - Reusable - Neonatal | | | • | • |
| 989803142711 | n/a | Verification gas | 1 gas cylinder | Mainstream CO ₂ sensor: M2501A, MP2: 865040; X2: 865039; X3: 867030 | • | • |

Capnography

Mainstream CO₂ solution



989803142651



989803142661



989803142671



989803206591



989803206591



989803142711

Cardiac output

Invasive cardiac output

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|--|--------|-----|
| 989803101031 | 23001A | Reusable injectate temp probe to measure cardiac output For use with Baxter CO-Set™ kit only Length: 2.4 m (7.9 ft) | 1 probe | C.O module: 862279; Hemodynamic ext.: 862111; Capnography ext.: 862187; Xper Flex Cardio: 860335, 860338 | | • |
| 989803101041 | 23001B | Reusable injectate temp probe to measure cardiac output For use with Baxter CO-Set™ kit only Length: 0.5 m (1.6 ft) | 1 probe | | | • |
| 989803101051 | 23002A | Ice bath temperature probe | 1 probe | CM12: 863303, 863366; CM120: 863302, 863321; CM150: 863304, 863322; XperFlex Cardio: 860335, 860338; Hemodynamic ext.: 867039, 862111; Capnography Ext.: 867040; Microstream ext.: 867041 | | • |
| 989803104611 | M1642A | Cardiac output cable for use with transpulmonary thermodilution Length: 2,9m (9.5 ft) | 1 cable | CM12: 863303, 863366; CM120: 863302; CM150: 863304; C.O. module: 862279; Hemodynamic ext.: 862111 | • | • |
| 989803104621 | M1643A | Cardiac output cable for use with pulse contour analysis and transpulmonary thermodilution Length: 4.8 m (15.7 ft) | 1 cable | C.O. module: 862279; Hemodynamic ext.: 862111; Capnography ext.: 862187; XPER Flex Cardio: 860335, 860338 | • | • |



989803101041
23001B

989803101051
23002A

989803104611
M1642A

989803104621
M1643A



Electrodes

Single-patient use adult ECG monitoring electrodes

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|---|---|--------|-----|
| 989803101671 | 40493D | Wet gel ECG electrode, adult, foam Shelf life: 6 months, diameter: 55 mm (2.13 in) | 5 electrodes/pouch 15 pouches/box 4 boxes/case (300 electrodes/case) | Adult ECG monitoring: IntelliVue, Avalon TOCO+ Transducer: 867249; Avalon Cableless Toco+ MP Transducer: 866075; Avalon Cableless ECG/IUP Transducer: 866077 | • | • |
| 989803101681 | 40493E | Wet gel ECG electrode, adult, foam Shelf life: 6 months, diameter: 54 mm (2.13 in) | 30 electrodes/pouch 10 pouches/box (300 electrodes/box) | | • | • |
| 989803105971 | M2202A | Radiotranslucent foam electrode, adult Shelf life: 6 months, diameter: 50 mm (1.97 in) | 30 electrodes/pouch 10 pouches/box (300 electrodes/box) | Adult ECG monitoring: IntelliVue | • | • |
| 989803129091 | M4612A | Solid gel ECG electrode, adult, foam, teardrop shape Shelf life: 6 months, diameter: 50 mm (1.97 in) | 5 electrodes/pouch 60 pouches/box (300 electrodes/box) | | • | • |
| 989803129101 | M4613A | Solid gel ECG electrode, adult, foam, teardrop shape Shelf life: 6 months, diameter: 50 mm (1.97 in) | 30 electrodes/pouch 10 pouches/box (300 electrodes/box) | | • | • |
| 989803148821 | | Solid gel ECG electrode, small adult, radiolucent foam Shelf life: 6 months, size: 35 x 40 mm (1.4 x 1.6 in) | 30 electrodes/pouch 20 pouches/box (600 electrodes/box) | Adult ECG monitoring: IntelliVue and Efficia | • | • |
| 989803192551 | | Wet gel ECG electrode, adult, clear tape | 30 electrodes/pouch 10 pouches/box (300 electrodes/box) | Adult ECG monitoring: IntelliVue | • | • |

Electrodes

Single-patient use adult ECG monitoring electrodes



Electrodes

Single-patient use adult metallic ECG lead wire electrodes – AAMI color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|--|----------------------------------|--------|-----|
| 989803156201 | | Pre-attached metallic 3-lead wire (0.91 m / 36 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in) | 3 electrodes/set 1 set/pouch 100 pouches/box (300 electrodes/box) | Adult ECG monitoring: IntelliVue | • | • |
| 989803156221 | | Pre-attached radiolucent 3-lead wire (0.91 m / 36 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in) | 3 electrodes/set 1 set/pouch 100 pouches/box (300 electrodes/box) | | • | • |
| 989803156241 | | Pre-attached metallic 5-lead wire (0.91 m / 36 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in) | 5 electrodes/set 1 set/pouch 60 pouches/box (300 electrodes/box) | | • | • |
| 989803156261 | | Pre-attached radiolucent 5-lead wire (0.91 m / 36 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in) | 5 electrodes/set 1 set/pouch 60 pouches/box (300 electrodes/box) | | • | • |



989803156201



989803156241


Electrodes

Single-patient use adult radiolucent ECG lead wire electrodes – IEC color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|--|----------------------------------|--------|-----|
| 989803156231 | | Pre-attached radiolucent 3-lead wire (1.0 m / 39.4 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in) | 3 electrodes/set 1 set/pouch 100 pouches/box (300 electrodes/box) | Adult ECG monitoring: IntelliVue | • | |
| 989803156271 | | Pre-attached radiolucent 5-lead wire (1.0 m / 39.4 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in) | 5 electrodes/set 1 set/pouch 60 pouches/box (300 electrodes/box) | | | • |



989803156231




989803156271

Electrodes

Single-patient use pediatric ECG monitoring electrodes

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|---|---|--------|-----|
| 989803100491 | 13951C | Disposable, repositionable solid gel electrode, snap 3 electrodes mounted on plastic card Shelf life: 6 months, dimensions: 22 x 33 mm (0.87 x 1.30 in) | 3 electrodes/pouch 100 pouches/box (300 electrodes/box) | Pediatric/neonatal ECG monitoring: IntelliVue | • | • |
| 989803100581 | 13955C | Disposable, repositionable solid gel electrode, snap 3 electrodes mounted on plastic card Shelf life: 6 months, dimensions: 22 x 22 mm (0.87 x 0.87 in) | 3 electrodes/pouch 100 pouches/box (300 electrodes/box) | | • | • |



989803100491
13951C

989803100581
13955C

Electrodes

Single-patient use infant/pediatric metallic lead wire electrodes – AAMI color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|---|--|--------|-----|
| 989803100501 | 13952A | Pre-attached metallic 46 cm (18 in) lead wire electrode, non-shielded, solid gel, AAMI Shelf life: 6 months, dimensions: 22 x 33 mm (0.87 x 1.30 in); socket size 1.5 mm (0.06 in) | 3 sets/pouch 100 pouches/box (300 sets/box) | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 3-ch ECG module: 989803110391, 862275; 3-ch ECG/respiration module: M1002A, 862276; VM3: 863077; VM4: 863063; VM6: 863064, 863065; VM8: 863066, 863068 | • | |
| 989803100541 | 13953A | Pre-attached metallic 46 cm (18 in) lead wire electrode, non-shielded, solid gel, AAMI Shelf life: 6 months, dimensions: 22 x 22 mm (0.87 x 0.87 in); socket size 1.5 mm (0.06 in) | 3 sets/pouch 100 pouches/box (300 sets/box) | | • | |
| 989803183081 | | Mini NeoLead, metallic, pre-attached 61 cm (24 in) lead wire pediatric/neonatal ECG electrodes, hydrogel/hydrocolloid adhesive, AAMI, standard DIN connector Diameter: 19 mm (0.75 in) | 3 sets/pouch 50 pouches/box (150 sets/box) | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 3-ch ECG module: 862275; 3-ch ECG/respiration module: 862276; VM6: 863064, 863065; VM8: 863066, 863068; XPER Flex Cardio: 860335, 860338; Requires ECG trunk cable: M1668A or 989803170181 | • | |
| 989803183101 | | NeoLead, metallic, pre-attached 61 cm (24 in) lead wire pediatric/neonatal ECG electrodes, hydrogel/hydrocolloid adhesive, AAMI, standard DIN connector Diameter: 25 mm (1 in) | 3 sets/pouch 50 pouches/box (150 sets/box) | | • | |



989803100501
13952A



989803100541
13953A



989803183081



989803183101


Electrodes

Single-patient use infant/pediatric radiolucent lead wire electrodes – AAMI color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|---|---|--------|-----|
| 989803100511 | 13952B | Pre-attached radiolucent 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 33 mm (0.87 x 1.30 in); socket size: 1.5 mm (0.06 in) | 3 sets/pouch 100 pouches/box (300 sets/box) | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 3-ch ECG module: 989803110391, 862275; 3-ch ECG/respiration module: M1002A, 862276; VM3: 863077; VM4: 863063; VM6: 863064, 863065; VM8: 863066, 863068 | • | • |
| 989803100551 | 13953B | Pre-attached radiolucent 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 22 mm (0.87 x 0.87 in); socket size: 1.5 mm (0.06 in) | 3 sets/pouch 100 pouches/box (300 sets/box) | | • | • |
| 989803183121 | | NeoLead, radiolucent pre-attached 61 cm (24 in) lead wire pediatric/neonatal ECG electrodes, AAMI, standard DIN connector Diameter: 25 mm (1 in) | 3 sets/pouch 50 pouches/box (150 sets/box) | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 3-ch ECG module: 862275; 3-ch ECG/respiration module: 862276; VM6: 863064, 863065; VM8: 863066, 863068; XPER Flex Cardio: 860335, 860338; Requires ECG trunk cable: M1668A or 989803170181 | • | • |
| 989803183141 | | Micro NeoLead, radiolucent pre-attached 61 cm (24 in) lead wire pediatric/neonatal ECG electrodes, AAMI, standard DIN connector Dimensions: 25 x 9.5 mm (1 x 0.37 in) | 3 sets/pouch 50 pouches/box (150 sets/box) | | • | • |



989803100511
13952B



989803100551
13953B



989803183121




989803183141

Electrodes

Single-patient use infant/pediatric metallic lead wire electrodes – IEC color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|---|----------|--------|-----|
| 989803100521 | 13952D | Pre-attached metallic 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 33 mm (0.87 x 1.30 in); socket size: 1.5 mm (0.06 in) | 3 sets/pouch 100 pouches/box (300 sets/box) | | • | |
| 989803100561 | 13953D | Pre-attached metallic 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 22 mm (0.87 x 0.87 in); socket size: 1.5 mm (0.06 in) | 3 sets/pouch 100 pouches/box (300 sets/box) | | • | |



989803100521
13952D




989803100561
13953D

Electrodes

Single-patient use infant/pediatric radiolucent lead wire electrodes – IEC color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe USA |
|---------------|---------------------------|---|---|---|------------|
| 989803100571 | 13953E | Pre-attached radiolucent 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 22 mm (0.87 x 0.87 in); socket size: 1.5 mm (0.06 in) | 3 sets/pouch 100 pouches/box (300 sets/box) | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 3-ch ECG module: 989803110391, 862275; 3-ch ECG/respiration module: M1002A, 862276; VM3: 863077; VM4: 863063; VM6: 863064, 863065; VM8: 863066, 863068 | • |



989803100571
13953E



Cables and leads


Lead sets – AAMI color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|--|--------|-----|
| 989803125841 | M1968A | ICU, grabber, 5 limb leads, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; | • | • |
| 989803125861 | M1973A | OR, grabber, 5 limb leads, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft) | 1 each | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304 | • | • |
| 989803125881 | M1976A | ICU, grabber, 5 chest leads, reusable Length: 1 m (3.3 ft) | 1 each | | • | • |
| 989803125901 | M1979A | OR, grabber, 5 chest leads, reusable Length: 1 m (3.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; CM12: 863303, 863366; CM120: 863302; CM150: 863304. 5+5 trunk cable M1949A, M1663A | • | • |
| 989803144841 | M1532A | ICU, grabber, 4 leads, length: 1 m (3.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 6+4 trunk cable M1665A | • | • |
| 989803144911 | M1602A | ICU, snap, 5 chest leads, reusable Length: 1 m (3.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; CM12: 863303, 863366; CM120: 863302; CM150: 863304. 5+5 trunk cable M1949A, M1663A | • | • |
| 989803144991 | M1644A | ICU, snap, 5 limb leads, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft) | 1 each | | • | • |
| 989803145091 | M1671A | ICU, grabber, 3 leads, length: 1 m (3.3 ft) | 1 each | | • | • |
| 989803145111 | M1673A | ICU, snap, 3 leads, length: 1 m (3.3 ft) | 1 each | | • | • |
| 989803145131 | M1675A | OR, grabber, 3 leads, length: 1 m (3.3 ft) | 1 each | | • | • |
| 989803145161 | M1680A | ICU, grabber, 6 leads, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 6-lead trunk cable M1667A or 6+4 trunk cable M1665A | • | • |
| 989803173121 | | ICU, grabber, 3 leads, length: 1 m (3.3 ft) | 20 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033 | • | • |
| 989803173131 | | ICU, grabber, 5 limb leads, single-patient use Length: 1 m (3.3 ft) | 20 each | | • | • |

Cables and leads

Lead sets – AAMI color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|--|--------|-----|
| 989803184931 | | Standard 10-lead patient cable, AHA Length: 1 m (3.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033 | • | • |
| 989803192141 | | ICU, grabber, 5 chest leads, single-patient use Length: 1 m (3.3 ft) | 20 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 5+5 trunk cable M1949A, M1663A | • | • |
| 989803197511 | | 6-lead set, single-patient use, bedside, AAMI Length: 1 m (3.3 ft) | 20 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033 | • | • |

989803125841
M1968A989803125861
M1973A989803125881
M1976A989803125901
M1979A989803144841
M1532A

Cables and leads

Lead sets – AAMI color coding



989803144911
M1602A



989803144991
M1644A




989803145091
M1671A



989803145111
M1673A



989803145131
M1675A



989803145161
M1680A



989803173121



989803173131



989803192141



989803197511



Cables and leads


Lead sets – IEC color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe USA |
|---------------|---------------------------|---|-----------|---|------------|
| 989803125851 | M1971A | ICU, grabber, 5 limb leads, reusable Length: limb 1.6 m (5.3 ft), chest 1 m (3.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; | • |
| 989803125871 | M1974A | OR, grabber, 5 limb leads, reusable Length: limb 1.6 m (5.3 ft), chest 1 m (3.3 ft) | 1 each | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304. 5-lead trunk cable M1668A, OR trunk cable 989803170181 or 5+5 trunk cable M1949A, M1663A | • |
| 989803125891 | M1978A | ICU, grabber, 5 chest leads, reusable Length: 1 m (3.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304. 5+5 trunk cable M1949A, M1663A | • |
| 989803144851 | M1533A | ICU, grabber, 4 chest leads, reusable Length: 1 m (3.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 6-lead trunk cable M1667A or 6+4 trunk cable M1665A | • |
| 989803145001 | M1645A | ICU, snap, 5 limb leads, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304. 5-lead trunk cable M1668A, OR trunk cable 989803170181 or 5+5 trunk cable M1949A, M1663A | • |
| 989803145101 | M1672A | ICU, grabber, 3 leads, reusable Length: 1 m (3.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; | • |
| 989803145121 | M1674A | ICU, snap, 3 leads, reusable Length: 1 m (3.3 ft) | 1 each | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304 | • |
| 989803145141 | M1678A | OR, grabber, 3 leads, reusable Length: 1 m (3.3 ft) | 1 each | | • |
| 989803145171 | M1681A | ICU, grabber, 6 leads, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft) | 1 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 6-lead trunk cable M1667A or 6+4 trunk cable M1665A | • |

Cables and leads

Lead sets – IEC color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe USA |
|---------------|---------------------------|---|-----------|---|------------|
| 989803174201 | | ICU, grabber, 3 leads, single-patient use Length: 1 m (3.3 ft) | 20 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304 | • |
| 989803174211 | | ICU, grabber, 5 limb leads, single-patient use Length: 1 m (3.3 ft) | 20 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304. 5-lead trunk cable M1668A, OR trunk cable 989803170181 or 5+5 trunk cable M1949A, M1663A | • |
| 989803192151 | | ICU, grabber, 5 chest leads, single-patient use Length: 1 m (3.3 ft) | 20 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 5+5 trunk cable M1949A, M1663A | • |
| 989803197521 | | ICU, grabber, 6 leads, single-patient use Length: 1 m (3.3 ft) | 20 each | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033 | • |



989803125851
M1971A



989803125871
M1974A



989803125891
M1978A



989803144851
M1533A



989803145001
M1645A

Cables and leads

Lead sets – IEC color coding



989803145101
M1672A



989803145121
M1674A



989803145141
M1678A



989803145171
M1681A



989803174201



989803174211



989803192151




989803197521

Cables and leads

IntelliVue cable family – leads accessories

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---------------------------------|--------------|------------|--------|-----|
| 989803148841 | | IntelliVue cable management kit | 5 organizers | All cables | • | • |



989803148841

Cables and leads

IntelliVue Telemetry System lead sets – AAMI color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|------------------------------------|--------|-----|
| 989803173151 | | Grabber and shower shield, 5 leads, single-patient use, length: 85 cm (34 in) | 20 each | TRx4851A: 862231; TRx4841A: 862439 | • | • |



989803173151

Cables and leads

IntelliVue Telemetry System shields and alignment guides

Miscellaneous MX40 cables/supplies

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|------------------------------|--------|-----|
| 989803172221 | | Patient monitor ECG trunk cable, AAMI/IEC, reusable, length: 2.7 m (9 ft) | 1 each | MX40: 865350, 865351, 865352 | • | • |



989803172221

Cables and leads

Efficia 3-lead sets

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|------------------------------|-----------|--|--------|-----|
| 989803160651 | | Efficia 3-lead grabber, AAMI | 1 each | CM10: 863301, 863365; CM12: 863303, 863366; | • | |
| 989803160661 | | Efficia 3-lead grabber, IEC | 1 each | CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322 | • | |

Efficia 5-lead sets

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|------------------------------------|-----------|--|--------|-----|
| 989803160691 | | Efficia 5-lead grabber, limb, AAMI | 1 each | CM10: 863301, 863365; CM12: 863303, 863366; | • | |
| 989803160701 | | Efficia 5-lead grabber, limb, IEC | 1 each | CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322 | • | |

Efficia 3/5-lead sets

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|---|--------|-----|
| 989803160641 | | Efficia 3/5 ECG Trunk Cable, AAMI/IEC, | 1 each | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322 | • | |



Cables and leads

Patient cables – AAMI color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|------------------------------|--------|-----|
| 989803171811 | | ECG + SpO ₂ , 3 leads, grabbers, reusable, length: 85 cm (34 in) | 1 each | MX40: 865350, 865351, 865352 | • | • |
| 989803171821 | | ECG, 5 leads, snaps, reusable, length: 85 cm (34 in) | 1 each | | • | • |
| 989803171831 | | ECG, 5 leads, grabbers, reusable, length: 85 cm (34 in) | 1 each | | • | • |
| 989803171841 | | ECG + SpO ₂ , 5 leads, snaps, reusable, length: 85 cm (34 in) | 1 each | | • | • |
| 989803171851 | | ECG + SpO ₂ , 5 leads, grabbers, reusable, length: 85 cm (34 in) | 1 each | | • | • |
| 989803171871 | | ECG + SpO ₂ , 6 leads, grabbers, reusable, length: 85 cm (34 in) | 1 each | | • | • |
| 989803172031 | | ECG, 5 leads, grabbers, single-patient use, length 85 cm (34 in) | 20 each | | • | • |
| 989803172051 | | ECG + SpO ₂ , 5 leads, grabbers, single-patient use, length: 85 cm (34 in) | 20 each | | • | • |
| 989803202411 | | ECG 5 leads, grabber AAMI, Masimo SpO ₂ , telemetry Length: 85 cm (34 in) | 1 each | | • | • |

Cables and leads

Patient cables – AAMI color coding



Cables and leads

Patient cables – IEC color coding

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|------------------------------|--------|-----|
| 989803171911 | | ECG + SpO ₂ , 3 leads, grabbers, reusable, length: 85 cm (34 in) | 1 each | MX40: 865350, 865351, 865352 | • | |
| 989803171931 | | ECG, 5 leads, grabbers, reusable, length: 85 cm (34 in) | 1 each | | • | |
| 989803171951 | | ECG + SpO ₂ , 5 leads, grabbers, reusable, length: 85 cm (34 in) | 1 each | | • | |
| 989803171971 | | ECG + SpO ₂ , grabbers, reusable, length: 85 cm (34 in) | 1 each | | • | |
| 989803172131 | | ECG, 5 leads, grabbers, single-patient use, length: 85 cm (34 in) | 20 each | | • | |
| 989803172151 | | ECG + SpO ₂ , 5 leads, grabbers, single-patient use, length: 85 cm (34 in) | 20 each | | • | |
| 989803202421 | | ECG 5 leads, grabber, IEC, Masimo SpO ₂ , telemetry length: 85 cm (34 in) | 1 each | | • | |

Cables and leads

Patient cables – IEC color coding



989803171911



989803171931



989803171951



989803171971



989803172131



989803172151

Cables and leads

Trunk cables

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|--|--------|-----|
| 989803125831 | M1949A | ECG trunk cable, 12-pin connector, AAMI and IEC labels included 10 (5+5) leads, reusable Length: 2.7 m (9.9 ft) | 1 cable | MP2: 865040; MP5: 865024; MMS: 862442; X2: 865039; X3: 867030; MX100: 867033; CM12: 863303, 863366; CM120: 863302, 863321; CM150: 863304, 863322; HeartStart XL: 861266; HeartStart XL+: 861290; HeartStart MRx: 861288, 861289; Xper Flex Cardio: 860338, 850335 | • | • |
| 989803144791 | M1663A | ECG trunk cable, 12-pin connector, AAMI and IEC labels included 10 (5+5) leads, reusable Length: 2 m (6.5 ft) | 1 cable | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 3 Ch ECG module: 862275; 3 ECG ch./respiration module: 862276; CM120: 863302, 863321; CM150: 863304, 863322; HeartStart XL: 861266; HeartStart XL+: 861290; HeartStart MRx: 861288 | • | • |
| 989803145041 | M1665A | ECG trunk cable, 12-pin connector, AAMI and IEC labels included 10 (6+4) leads, reusable Length: 2.7 m (9.9 ft) | 1 cable | MP2: 865040; MP5: 865024; MMS: 862442; X2: 865039; X3: 867030; MX100: 867033 | • | • |
| 989803145051 | M1667A | ECG trunk cable, 12-pin connector, AAMI and IEC labels included 6 leads, reusable Length: 2.7 m (9.9 ft) | 1 cable | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; HeartStart XL+: 861290 | • | • |
| 989803145061 | M1668A | ECG trunk cable, 12-pin connector, AAMI and IEC labels included 5 leads, reusable Length: 2.7 m (9.9 ft) | 1 cable | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 3-ch ECG module: 862275; 3-ch ECG/respiration module: 862276; VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088, 863317; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322; DFM100: 866199; HeartStart XL: 861266; HeartStart XL+: 861290; HeartStart MRx: 861288, 861289; Xper Flex Cardio: 860338, 850335 | • | • |

Cables and leads

Trunk cables

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|--|--------|-----|
| 989803145071 | M1669A | ECG trunk cable, 12-pin connector, AAMI and IEC labels included 3 leads, reusable Length: 2.7 m (9.9 ft) | 1 cable | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 3-ch ECG module: 862275; 3-ch ECG/respiration module: 862276; VM3: 863077; VM4: 863063, 863085; VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088, 863317; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322; DFM100: 866199; HeartStart XL: 861266; HeartStart XL+: 861290; HeartStart MRx: 861288, 861289 | • | • |
| 989803170171 | | ECG trunk cable, orange, 12-pin connector, AAMI and IEC labels included For use in operating room 3 leads, reusable Length: 2.7 m (9.9 ft) | 1 cable | MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; VM6: 863064, 863065; VM8: 863066, 863068, 863317; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199 | • | • |
| 989803170181 | | ECG trunk cable, orange, 12-pin connector, AAMI and IEC labels included For use in operating room 5 leads, reusable Length: 2.7 m (9.9 ft) | 1 cable | MP2: 865040; MP5: 865024; MMS: 862442; X2: 865039; X3: 867030; MX100: 867033; VM6: 863065; VM8: 863066, 863068, 863317; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322; DFM100: 866199 | • | • |

ECG

Cables and leads

Trunk cables



989803125831
M1949A



989803144791
M1663A



989803145041
M1665A



989803145051
M1667A



989803145061
M1668A



989803145071
M1669A



989803170171



989803170181



Telemetry adapters

ECG telemetry adapters

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|------------------------------|--------|-----|
| 989803199071 | | MX40 telemetry adapter, ECG + SpO ₂ , 6 leads | 1 each | MX40: 865350, 865351, 865352 | • | • |
| 989803199081 | | MX40 telemetry adapter, SpO ₂ only | 1 each | | • | • |
| 989803199091 | | MX40 telemetry adapter, ECG + SpO ₂ , 5 leads | 1 each | | • | • |
| 989803199101 | | MX40 telemetry adapter, ECG + SpO ₂ , 3 leads | 1 each | | • | • |
| 989803202381 | | MX40 Adapter 5-lead + Masimo SpO ₂ | 1 each | | • | • |
| 989803202391 | | MX40 Adapter Masimo SpO ₂ -only | 1 each | | • | • |
| 989803204291 | | MX40 Adapter 6-lead + Masimo SpO ₂ | 1 each | | • | • |



989803199071



989803199081



989803199091



989803199101



989803202381




989803202391

EEG supplies

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|----------------------------|--------|-----|
| 989803105621 | M1931A | 5-lead EEG adult cup electrode set, reusable, length: 80 cm (31.5 in) | 5 leads | EEG module: 862290, 866078 | • | • |
| 989803105631 | M1932A | 5-lead EEG pediatric cup electrode set, reusable, length: 80 cm (31.5 in) | 5 leads | | • | • |
| 989803105641 | M1934A | 5-lead EEG miniclip lead set, reusable, length: 80 cm (31.5 in) | 5 leads | | • | • |
| 989803106091 | M2268A | 5-lead trunk cable, reusable, length: 2.7 m (8.8 ft) | 1 cable | EEG module: 862290, 866078 | • | • |
| 989803180521 | | 9-lead EEG infant and pediatric cup electrode set, reusable, length: 80 cm (31.5 in) | 9 leads | | • | • |
| 989803180531 | | 9-lead EEG miniclip lead set, reusable, length: 80 cm (31.5 in) | 9 leads | | • | • |
| 989803180541 | | 4-channel trunk cable, reusable, length: 2.7 m (8.8 ft) | 1 cable | | • | • |
| 989803180561 | | 9-lead EEG adult cup electrode set, reusable, length: 80 cm (31.5 in) | 9 leads | | • | • |

EEG supplies



989803105621
M1931A



989803105631
M1932A



989803105641
M1934A



989803106091
M2268A



989803180521



989803180531



989803180541




989803180561

Electrodes

Single-patient use pediatric EEG monitoring electrodes

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|---|--------------------|--------|-----|
| 989803100491 | 13951C | Disposable, repositionable solid gel electrode, snap 3 electrodes mounted on plastic card Shelf life: 6 months, dimensions: 22 x 33 mm (0.87 x 1.30 in) | 3 electrodes/pouch 100 pouches/box (300 electrodes/box) | EEG module: 862290 | • | • |
| 989803100581 | 13955C | Disposable, repositionable solid gel electrode, snap 3 electrodes mounted on plastic card Shelf life: 6 months, dimensions: 22 x 22 mm (0.87 x 0.87 in) | 3 electrodes/pouch 100 pouches/box (300 electrodes/box) | | • | • |



989803100491
13951C



989803100581
13955C

NMT

IntelliVue NMT patient cable

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|---------------|--------------------|--------|-----|
| 989803174581 | | Patient cable with 12-pin instrument connector, reusable, length: 3.0 m (9.8 ft) | 1 each | NMT module: 865383 | • | • |
| 989803199211 | | IntelliVue NMT hand adapter | 10 pieces/box | | • | • |



989803174581



989803199211

Non-invasive blood pressure

Multi-patient cuffs

Comfort Care cuffs – single hose

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|--|--------|-----|
| 989803104141 | M1571A | Comfort Care cuff, infant Bladder length: 15 cm (5.9 in) Bladder width: 5.5 cm (2.2 in) Limb circumference: 10–15 cm (3.9–5.9 in) | 1 cuff | NBP module: 989803110421, 862278; Series 5000: 862318; Series 50 FM: 862334, 862335; FM20: 862198 (except M4572B/4573B); FM30: 862199 (except M4572B/4573B); FM40: 865132 (except M4572B/4573B); FM50: 865071 (except M4572B/4573B); Xper Flex Cardio: 860335, 860338; VM4: 863063; VM6: 863064, 863065; VM8: 863066, 863068; VS3: 863069, 863070, 863071, 863072, 863073, 863074; VS2: 863079, 863080, 863081, 863082, 863278, 863279; VSi: 863275, 863276, 863277; VS4: 863283; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; MP2: 865040; MP5: 865024, 865322, 865120; MMS: 862442; X2: 865039; X3: 867030; MX100: 867033; HeartStart XL+: 861290; HeartStart MRx: 861288, 861493, 861494, 861495, 861289, 861464, 861465, 861480, 861481, 861482, 861483, 861484, 861491, 860396, 860397, 860398 | • | • |
| 989803104151 | M1572A | Comfort Care cuff, pediatric Bladder length: 21.5 cm (8.5 in) Bladder width: 8 cm (3.1 in) Limb circumference: 14–21.5 cm (5.5–8.5 in) | 1 cuff | | • | • |
| 989803104161 | M1573A | Comfort Care cuff, small adult Bladder length: 24 cm (9.4 in) Bladder width: 10.5 cm (4.1 in) Limb circumference: 20.5–28 cm (8–11 in) | 1 cuff | | • | • |
| 989803104171 | M1574A | Comfort Care cuff, adult Bladder length: 28 cm (11 in) Bladder width: 13 cm (5.1 in) Limb circumference: 27–35 cm (10.6–13.7 in) | 1 cuff | | • | • |
| 989803104181 | M1575A | Comfort Care cuff, large adult Bladder length: 35 cm (13.8 in) Bladder width: 16 cm (6.3 in) Limb circumference: 34–43 cm (13.3–16.9 in) | 1 cuff | | • | • |



989803104141
M1571A



989803104151
M1572A



989803104161
M1573A



989803104171
M1574A



989803104181
M1575A

Non-invasive blood pressure

Multi-patient cuffs

Comfort Care cuffs – single hose

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|---|--------|-----|
| 989803104191 | M1576A | Comfort Care cuff, large adult XL Bladder length: 43 cm (16.9 in) Bladder width: 20 cm (7.9 in) Limb circumference: 42–54 cm (16.5–21.2 in) | 1 cuff | NBP module: 989803110421, 862278; Series 5000: 862318; Series 50 FM: 862334, 862335; FM20: 862198; FM30: 862199; FM40: 865132; FM50: 865071; Xper Flex Cardio: 860335, 860338; VM4: 863063; VM6: 863064, 863065; VM8: 863066, 863068; VS3: 863069, 863070, 863071, 863072, 863073, 863074; VS2: 863079, 863080, 863081, 863082, 863278, 863279; VSi: 863275, 863276, 863277; VS4: 863283; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; MP2: 865040; MP5: 865024, 865322, 865120; MMS: 862442; X2: 865039; X3: 867030; MX100: 867033; HeartStart XL+: 861290; HeartStart MRx: 861288, 861493, 861494, 861495, 861289, 861464, 861465, 861480, 861481, 861482, 861483, 861484, 861491, 860396, 860397, 860398 | • | • |
| 989803191311 | M1574XL | Comfort Care cuff, adult XL Bladder length: 28 cm (11 in) Bladder width: 13 cm (5.1 in) Limb circumference: 27–35 cm (10.6–13.7 in) | 1 cuff | | • | • |



989803104191
M1576A

989803191311
M1574XL

Non-invasive blood pressure

Multi-patient cuffs

Easy Care cuffs – single hose

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|----------|--------|-----|
| 989803147841 | M4552B | Infant, orange Bladder length: 15 cm (5.9 in) Bladder width: 5.5 cm (2.2 in) Limb circumference: 10–15 cm (3.9–5.9 in) | 1 cuff | | • | • |
| 989803147851 | M4553B | Pediatric, green Bladder length: 21.5 cm (8.5 in) Bladder width: 8 cm (3.1 in) Limb circumference: 14–21.5 cm (5.5–8.5 in) | 1 cuff | | • | • |
| 989803147861 | M4554B | Small adult, light blue Bladder length: 24 cm (9.4 in) Bladder width: 10.6 cm (4.2 in) Limb circumference: 20.5–28.5 cm (8–11.2 in) | 1 cuff | | • | • |
| 989803147871 | M4555B | Adult, navy blue Bladder length: 30 cm (11.8 in) Bladder width: 13.5 cm (5.3 in) Limb circumference: 27.5–36 cm (10.8–14.1 in) | 1 cuff | | • | • |

Non-invasive blood pressure

Multi-patient cuffs


Easy Care cuffs – single hose

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|--|--------|-----|
| 989803147881 | M4556B | Adult extra-long, navy blue Bladder length: 30 cm (11.8 in) Bladder width: 13.5 cm (5.3 in) Limb circumference: 27.5–36 cm (10.8–14.1 in) | 1 cuff | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; NBP module: 862278; MMS: 862392, 862442; X2: 865039; MP2: 865040; MP5: 865024, 865120, 865322; MB1 bundle: 866044; MB2 bundle: 866045; 865322; MB1 bundle: 866044; MB2 bundle: 866045; | • | • |
| 989803147891 | M4557B | Large adult, burgundy Bladder length: 38 cm (15 in) Bladder width: 17 cm (6.7 in) Limb circumference: 35–45 cm (13.7–17.7 in) | 1 cuff | 865322; MB1 bundle: 866044; MB2 bundle: 866045; MB3 bundle: 866046; C3: 862474, 862478; VM4: 863063, 863085, 863287; VM6: 863064, 863065, 863086, 863288; VM8: 863066, 863068, 863087, 863088, 863289; VS1: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863062; VS2: 863079, 863080, 863081, 863082, 863278, 863279; VS3: 863069, 863070, 863071, 863072, 863073, 863074, 863285; VS4: 863283, 863286; VS1: 863275, 863276, 863277; FM20: 862198; FM30: 862199; FM40: 865132; FM50: 865071; Series 50: 852334, 852335, 852336, 852337 | • | • |
| 989803147901 | M4558B | Large adult extra-long, burgundy Bladder length: 38 cm (15 in) Bladder width: 17 cm (6.7 in) Limb circumference: 35–45 cm (13.7–17.7 in) | 1 cuff | | • | • |
| 989803147911 | M4559B | Thigh, grey Bladder length: 48 cm (18.9 in) Bladder width: 21 cm (8.3 in) Limb circumference: 44–56 cm (17.3–22 in) | 1 cuff | | • | • |

Non-invasive blood pressure

Multi-patient cuffs

Easy Care cuffs – single hose





Non-invasive blood pressure

Single-patient cuffs Single Care cuffs

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe USA |
|---------------|---------------------------|--|-----------|---|------------|
| 989803182281 | | Pediatric, green | 20 cuffs | NBP module: 862278; MMS: 862442, 866462; | • |
| 989803182291 | | Small adult, light blue | 20 cuffs | MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MB1 bundle: 866044; MB2 bundle: 866045; MB3 bundle: 866046; | • |
| 989803182301 | | Adult, navy blue | 20 cuffs | VS1: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863062; C3: 862474, 862478; VM4: 863063, 863085, 863287; VM6: 863064, 863065, 863086, 863288; VM8: 863066, 863068, 863087, 863088, 863289; VS2: 863079, 863080, 863081, 863082, 863278, 863279; VS3: 863069, 863070, 863071, 863072, 863073, 863074, 863285; VS4: 863283, 863286; VSi: 863275, 863276, 863277; | • |
| 989803182311 | | Adult XL, navy blue with edge striping | 20 cuffs | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; XPER Flex Cardio: 860335, 860338 | • |
| 989803182321 | | Large adult, burgundy | 20 cuffs | | |

Available for 989803182291 / 989803182301 / 989803182311 / 989803182321:
FM20: 862198; FM30: 862199; FM40: 865132; FM50: 865071



989803182281

989803182291

989803182301

989803182311

989803182321

Non-invasive blood pressure

Single-patient cuffs

Soft neonatal/infant cuffs – single hose

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|---|---|--------|-----|
| 989803167231 | M1866S | Neonatal soft single-patient cuff, size #1, single hose | 40 cuffs | NBP module: 862278; MMS: 862442, 866462; MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MB1 bundle: 866044; MB2 bundle: 866045; MB3 bundle: 866046; VS1: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863062; C3: 862474, 862478; VM4: 863063, 863085, 863287; VM6: 863064, 863065, 863086, 863288; VM8: 863066, 863068, 863087, 863088, 863289; VS2: 863079, 863080, 863081, 863082, 863278, 863279; VS3: 863069, 863070, 863071, 863072, 863073, 863074, 863285; VS4: 863283, 863286; VS5: 863275, 863276, 863277; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304 | • | • |
| 989803167241 | M1868S | Neonatal soft single-patient cuff, size #2, single hose | 40 cuffs | | • | • |
| 989803167251 | M1870S | Neonatal soft single-patient cuff, size #3, single hose | 40 cuffs | | • | • |
| 989803167261 | M1872S | Neonatal soft single-patient cuff, size #4, single hose | 40 cuffs | | • | • |
| 989803167271 | M1873S | Infant soft single-patient cuff, size #5, single hose | 40 cuffs | | • | • |
| 989803167551 | | Neonatal soft single-patient cuffs, single hose | 50 cuffs Consists of: 5 x M1866S, 10 x M1868S, 20 x M1870S, 10 x M1872S & 5 x M1873S | | • | • |

Non-invasive blood pressure

Single-patient cuffs

Soft neonatal/infant cuffs – single hose



989803167231
M1866S



989803167241
M1868S



989803167251
M1870S



989803167261
M1872S



989803167271
M1873S



989803167551

Non-invasive blood pressure

Single-patient cuffs

Vinyl neonatal/infant cuffs

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|---|---|--------|-----|
| 989803167181 | M1866B | Neonatal single-patient cuff, size 1, single hose New connector configuration, replaces the M1866A Limb circumference: 3.1–5.7 cm (1.2–2.2 in) | 40 cuffs | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; NBP module: 862278; MMS: 862392, 862442; X2: 865039; MP2: 865040; MP5: 865024, 865322; VM4: 863063; VM6: 863064, 863065; VM8: 863066, 863068; VS1: 863055, 863056, 863057, 863058, 863059, 863060, 863061; VS2: 863079, 863080, 863081, 863082, 863278; VS3: 863069, 863070, 863071, 863072, 863073, 863074; VSi: 863275, 863276; Legacy HP monitors: 78352A, 78352C, 78834C | • | • |
| 989803167191 | M1868B | Neonatal single-patient cuff, size 2, single hose New connector configuration, replaces the M1868A Limb circumference: 4.3–8.0 cm (1.7–3.1 in) | 40 cuffs | | • | • |
| 989803167201 | M1870B | Neonatal single-patient cuff, size 3, single hose New connector configuration, replaces the M1870A Limb circumference: 5.8–10.9 cm (2.2–4.3 in) | 40 cuffs | | • | • |
| 989803167211 | M1872B | Neonatal single-patient cuff, size 4, single hose New connector configuration, replaces the M1872A Limb circumference: 7.1–13.1 cm (2.8–5.1 in) | 40 cuffs | | • | • |
| 989803167221 | M1873B | Infant single-patient cuff, size 5, single hose New connector configuration Limb circumference: 10–15 cm (3.9–5.9 in) | 40 cuffs | | • | • |
| 989803167541 | | Neonatal single-patient NBP kit New connector configuration, replaces the M1820-60020 | 50 cuffs Consists of: 5x M1866B size 1, 10x M1868B size 2, 20x M1870B size 3, 10x M1872B size 4, 5x M1873B size 5 | • | • | |

Non-invasive blood pressure

Single-patient cuffs

Vinyl neonatal/infant cuffs



989803167181
M1866B
989803167191
M1868B
989803167201
M1870B
989803167211
M1872B
989803167221
M1873B

989803167541

Non-invasive blood pressure

Single-patient cuffs Gentle Care cuffs – single hose

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|---|--------|-----|
| 989803148001 | M4572B | Infant, orange Bladder length: 15 cm (5.9 in) Bladder width: 5.5 cm (2.2 in) Limb circumference: 10–15 cm (3.9–5.9 in) | 10 cuffs | NBP module: 862278; MMS: 862442, 866462; MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MB1 bundle: 866044; MB2 bundle: 866045; MB3 bundle: 866046; VS1: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863062; C3: 862474, 862478; VM4: 863063, 863085, 863287; VM6: 863064, 863065, 863086, 863288; VM8: 863066, 863068, 863087, 863088, 863289; VS2: 863079, 863080, 863081, 863082, 863278, 863279; VS3: 863069, 863070, 863071, 863072, 863073, 863074, 863285; VS4: 863283, 863286; VSi: 863275, 863276, 863277; | • | • |
| 989803148011 | M4573B | Pediatric, green Bladder length: 21.5 cm (8.5 in) Bladder width: 8 cm (3.1 in) Limb circumference: 14–21.5 cm (5.5–8.5 in) | 10 cuffs | | • | • |
| 989803148021 | M4574B | Small adult, light blue Bladder length: 24 cm (9.4 in) Bladder width: 10.6 cm (4.2 in) Limb circumference: 20.5–28.5 cm (8–11.2 in) | 10 cuffs | | • | • |
| 989803148031 | M4575B | Adult, navy blue Bladder length: 30 cm (11.8 in) Bladder width: 13.5 cm (5.3 in) Limb circumference: 27.5–36 cm (10.8–14.1 in) | 10 cuffs | XPER Flex Cardio: 860335, 860338; FM20: 862198 (except M4572B/4573B); FM30: 862199 (except M4572B/4573B); FM40: 865132 (except M4572B/4573B); FM50: 865071 (except M4572B/4573B) | • | • |
| 989803148041 | M4576B | Adult extra-long, navy blue Bladder length: 30 cm (11.8 in) Bladder width: 13.5 cm (5.3 in) Limb circumference: 27.5–36 cm (10.8–14.1 in) | 10 cuffs | | • | • |
| 989803148051 | M4577B | Large adult, burgundy Bladder length: 38 cm (15 in) Bladder width: 17 cm (6.7 in) Limb circumference: 35–45 cm (13.7–17.7 in) | 10 cuffs | | • | • |
| 989803148061 | M4578B | Large adult extra-long, burgundy Bladder length: 38 cm (15 in) Bladder width: 17 cm (6.7 in) Limb circumference: 35–45 cm (13.7–17.7 in) | 10 cuffs | | • | • |

Non-invasive blood pressure

Single-patient cuffs Gentle Care cuffs – single hose



989803148001
M4572B
989803148011
M4573B
989803148021
M4574B
989803148031
M4575B
989803148041
M4576B
989803148051
M4577B
989803148061
M4578B



Non-invasive blood pressure

Accessories and hoses

Blood pressure air hoses

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|------------|---|--------|-----|
| 989803104331 | M1598B | Adult NBP air hose Length: 1.5 m (4.92 ft) | 1 air hose | NBP module: 862278; MMS: 862442, 866462; MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MB1 bundle: 866044; MB2 bundle: 866045; MB3 bundle 866046; VS1: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863062; C3: 862474, 862478; VM4: 863063, 863085, 863287; VM6: 863064, 863065, 863086, 863288; VM8: 863066, 863068, 863087, 863088, 863289; VS2: 863079, 863080, 863081, 863082, 863278, 863279; VS3: 863069, 863070, 863071, 863072, 863073, 863074, 863285; VS4: 863283, 863286; VSi: 863275, 863276, 863277; Series 50 FM: 862334, 862335, 862336, 862337; FM20: 862198; FM30: 862199; FM40: 865132; FM50: 865071; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304 Except for product M1598B: Xper Flex Cardio; 860335, 860338 | • | • |
| 989803104341 | M1599B | Adult NBP air hose Length: 3.0 m (9.8 ft) | 1 air hose | | • | • |
| 989803166851 | M1596C | Neonatal NBP air hose Length: 1.5 m (4.92 ft) | 1 air hose | | • | • |
| 989803166861 | M1597C | Neonatal NBP air hose Length: 3.0 m (9.8 ft) | 1 air hose | | • | • |



989803104331
M1598B



989803104341
M1599B



989803166851
M1596C



989803166861
M1597C

Invasive blood pressure

Reusable invasive blood pressure solution

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|--------------|--|--------|-----|
| 989803102481 | CPJ84022 | Single-use sterile domes | 50 domes | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; CPJ840J5; CPJ840J6; IBP module: M1006A, M1006B, MMS: 862392, 862442; X2: 865039; MP2: 865040; MP5: 865024; Hemodynamic ext.: 862111; Capnography ext.: 862187; Mainstream ext.: 862394; Microstream® ext.: 862393, 865377; VM6: 863065; VM8: 863066, 863068 | • | • |
| 989803102491 | CPJ84046 | Transducer holder | 4 holders | MP2: 865040; MP5: 865024, 865120; X2: 865039; X3: 867030; MX100: 867033; Microstream® ext.: 862393, 865377; HeartStart MRx: 861288, 861289, 861464, 861480, 861481, 861482, 861483, 861484 | • | • |
| 989803102511 | CPJ840J6 | Reusable pressure transducer, used to measure arterial and venous blood pressure Includes one transducer holder Cable length: 3 m (9.8 ft) | 1 transducer | MP2: 865040; MP5: 865024; IBP module: M1006A, M1006B; MMS: 862392, 862442; X2: 865039; X3: 867030; Capnography ext.: 862187; Mainstream CO ₂ module 862282; Mainstream ext.: 862394; VM6: 863065; VM8: 863066, 863068; Microstream® ext.: 862393, 865377; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304 | • | • |
| 989803102521 | CPJ84447 | IV pole mount | 1 mount | CPJ840J5 or CPJ840J6 with instruments: M8102A, M8105A; IBP module: M1006A, M1006B, M1312A, MMS: 862392, 862442; X2: 865039; MP2: 865040; MP5: 865024; Capnography ext.: 862187; Mainstream ext.: 862394; Microstream® ext.: 862393, 865377; VM6: 863065; VM8: 863066, 863068 | • | • |
| 989803136281 | CPJ84026 | Monitoring line set | 10 sets/bag | VM6: 863065; VM8: 863066, 863068, 863087, 863088; HeartStart MRx: 861288, 861289; X2: 865039; X3: 867030; MP2: 865040; MP5: 865024, 865120; Microstream® ext.: 862393, 865377 | • | • |
| 989803199741 | | Dual IBP adapter | 1 adapter | MX100: 867033; X3: 867030; Hemodynamic ext.: 867039; Capnography ext.: 867040; Microstream® ext.: 867041 | • | • |

Invasive blood pressure

Reusable invasive blood pressure solution



989803102481
CPJ84022



989803102491
CPJ84046



989803102511
CPJ840J6



989803102521
CPJ84447



989803136281
CPJ84026



989803199741

Paper

Recording paper

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|---|--------|-----|
| 989803101571 | 40477A | Thermal array recorder paper, 50 mm x 30.5 m (1.9 in x 100 ft) | 10 rolls | Recorder: M1116A, M1116B, M3176B; MP 20 Junior, MP20, MP30, MP40, MP50, MP60, MP70, MP90; MMS: 862442 | • | • |
| 989803109871 | M4816A | Thermal array recorder paper, 50 mm x 30.5 m (1.9 in x 100 ft) | 10 rolls | M3925A; M3176C; M1116C | • | • |
| 989803136891 | | Thermal array recorder paper, 5.8 cm x 30.5 m (2.2 in x 100 ft) | 5 rolls | VS3: 863069; VM4: 863063, 863085; VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088, 863317; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304 | • | • |
| 989803190381 | | Efficia chemical thermal recorder paper, 50 mm (1.9 in), no lines | 12 rolls | DFM 100: 866199 | • | |



989803101571
40477A



989803109871
M4816A



989803136891



989803190381

Pulse oximetry

Reusable sensors

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|---|--------|-----|
| 989803144371 | M1191B | Reusable adult SpO ₂ sensor, 8-pin connector Length: 2 m (6.5 ft) Patient size: > 50 kg (> 110 lb) | 1 sensor | Pulse oximetry module: 862112; MMS: 862442, 866462; MP2: 865040; MP5: 865024, 865120; X2: 865039; | • | • |
| 989803144381 | M1191BL | Reusable adult SpO ₂ sensor, 8-pin connector Length: 3 m (9.8 ft) Patient size: > 50 kg (> 110 lb) | 1 sensor | X3: 867030; MX40: 865350, 865351, 867146; MX100: 867033; MMX: 867036; VM1: 863264, 863265, 863266; VM4: 863063, 863085; VM6: 863064, 863065, 863086; | • | • |
| 989803128591 | M1191T | Reusable adult SpO ₂ sensor, 9-pin (D-sub) connector Length: 0.45 m (1.5 ft) Patient size: > 50 kg (> 110 lb) | 1 sensor | VM8: 863066, 863068, 863087, 863088, 863317; VS2: 863278, 863279; VS4: 863283; VSi: 863276; VS30: 863359; CM10: 863301, 863365; CM12: 863303, 863366; | • | • |
| 989803205871 | M1192A | Reusable pediatric/small adult finger SpO ₂ Sensor, 8-pin connector Length: 1.5 m (4.9 ft) Patient size: 15–50 kg (33–110 lb) | 1 sensor | CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322; | • | • |
| 989803205901 | M1196A | Reusable adult SpO ₂ clip sensor, 8-pin connector Length: 3 m (9.8 ft) Patient size: > 40 kg (> 88 lb) | 1 sensor | HeartStart MRx: 861288, 861289, 861464, 861482, 861491, 860396, 860397, 860398, 861465; HeartStart XL+: | • | • |
| 989803205911 | M1196T | Reusable adult SpO ₂ clip sensor, 9-pin (D-sub) connector Length: 0.90 m (2.9 ft) Patient size: >40 kg (> 88 lb) | 1 sensor | FM30: 862199; FM40: 865132; FM50: 865071. | • | • |
| | | | | Except product M1191T, M1196T: DFM100: 866199 | • | • |

Pulse oximetry

Reusable sensors

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|--|--------|-----|
| 989803205881 | M1193A | Reusable neonatal hand/foot SpO ₂ sensor, 8-pin connector Length: 1.5 m (4.9 ft) Patient size: 1–4 kg (2.2–8.8 lb) | 1 sensor | Pulse oximetry module: 862112; MMS: 862442, 866462; MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MMX: 867036; VM1: 863264, 863265, 863266; VM3: 863077; VM4: 863063, 863085; VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088, 863317; VS2: 863278, 863279; VS3: 863071, 863072, 863073, 863074; VS4: 863283; VSi: 863276; VS30: 863359; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322 HeartStart MRx: 861288, 861289, 861464, 861482, 861491, 860396, 860397, 860398, 861465 Except product M1193A: FM30: 862199; FM40: 865132; FM50: 865071; MX40: 865350, 865351, 867146; HeartStart XL+: 861290. | • | • |
| 989803205891 | M1194A | Reusable adult/pediatric ear clip SpO ₂ sensor, 8-pin connector Length: 1.5 m (4.9 ft) Patient size: > 40 kg (> 88 lb) | 1 sensor | | • | • |

Pulse oximetry

Reusable sensors

| | | | |
|--|--|--|--|
|  |  |  |  |
| 989803144371 M1191B | 989803144381 M1191BL | 989803128591 M1191T | 989803205871 M1192A |
|  |  |  |  |
| 989803205901 M1196A | 989803205911 M1196T | 989803205881 M1193A | 989803205891 M1194A |

Pulse oximetry

Single-patient sensors

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|------------|---|--------|-----|
| 989803205381 | | Nasal Alar SpO ₂ sensor, compatible with FAST technology | 20 sensors | Pulse oximetry module: 862112; MMS: 862442; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX40: 865350, 865351, 865352; MX100: 867033; MMX: 867036; MP2: 865040 FM30: 862199; FM40: 865132; FM50: 865071 | • | • |
| 989803205391 | | Nasal Alar SpO ₂ sensor, compatible with FAST and Nellcor™ technologies | 24 sensors | Pulse oximetry module: 862112; MMS: 862442; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX40: 865350, 865351, 865352; MX100: 867033; MMX: 867036; MP2: 865040 Nellcor N-600x; Nellcor N-595; DataScope Passport 2 FM30: 862199; FM40: 865132; FM50: 865071 | • | • |
| 989803205831 | M1131A | Disposable adult/pediatric SpO ₂ sensor, 9-pin (D-sub) connector Length: 45 cm (17.7 in) Patient size: > 20 kg (> 44 lb) | 20 sensors | Pulse oximetry module: 862112; MMS: 862442, MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX40: 865350, 865351, 867146; MX100: 867033; MMX: 867036; VM1: 863264, 863265, 863266; VM4: 863063, 863085; VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088, 863317; VS2: 863278, 863279; VSi: 863276, 863277; VS30: 863359 CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322 HeartStart MRx: 861288, 861289, 861464, 861482, 861491, 860396, 860397, 860398, 861465; HeartStart XL+: 861290; FM30: 862199; FM40: 865132; FM50: 865071; XPER: 860338; EarlyVue: 863380; Intrepid: 867172 | • | • |

Pulse oximetry

Single-patient sensors

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|------------|--|--------|-----|
| 989803205851 | M1133A | Neonatal/infant/adult disposable SpO ₂ sensor, 9-pin (D-sub) connector Length: 0.90 m (2.9 ft) Patient sizes: Adult (finger): > 40 kg (> 88 lb) Neonate (foot/hand): < 3 kg (< 6.6 lb) Infant (big toe/thumb): 10–20 kg (22–44 lb) | 20 sensors | Pulse oximetry module: 862112; MMS: 862442; MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MMX: 867036; VM1: 863264, 863265, 863266; VM4: 863063, 863085; VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088, 863317; VS2: 863278, 863279; VS4: 863283; VSi: 863276; 863277; VS30: 863359, 863380; Intrepid: 867172; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322, M1132A | • | • |
| 989803205861 | M1134A | Adhesive-free neonatal/infant/adult disposable SpO ₂ sensor, 9-pin (D-sub) connector Length: 0.90 m (2.9 ft) Patient sizes: Adult (finger): > 40 kg (> 88 lb) Neonate (foot/hand): < 3 kg (< 6.6 lb) Infant (big toe/thumb): 10–20 kg (22–44 lb) | 20 sensors | Except product M1134A: 860338; HeartStart XL+: 861290. | • | • |
| 989803205841 | M1132A | Infant disposable SpO ₂ sensor, 9-pin (D-sub) connector Length: 0.9 m (2.9 ft) Patient sizes: Infant (finger/toe): 3-10kg (6.6-22 lb) | 20 sensors | Except M1133A/M1134A: 860335; FM30: 862199; FM40: 865132; FM50: 865071; MX40: 865350, 865351, 867146. | | |



989803205381; 989803205391
989803205401



989803205831
M1131A



989803205841
M1132A



989803205861; 989803205851
M1133A; M1134A

Adapter cables

Adapter cables for Philips sensors

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|---|--------|-----|
| 989803105681 | M1941A | SpO ₂ Extension Cable (8-pin), adds additional length between 8-pin sensors and 8-pin sockets Length: 2 m (6.6 ft) | 1 cable | Pulse oximetry module: 862112; MMS: 862442; MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MMX: 867036; VM4: 863063, 863085; VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088, 863317; VS2: 863278, 863279; VS4: 863283; VSi: 863276; VS30: 863359; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322; HeartStart XL+: 861290 FM30: 862199; FM40: 865132; FM50: 865071 Except product M1941A, M1943A/AL: 863277, 865350, 865351, 867146, 867172, M2703A, M2704A, M2705A Except product M1941A only: 863310, 863311, 863276, 863359 Except product M1943A only: 860335, 860338 Except product M1943A/AL only: 863380 Except product M1943A, M1943AL, M1143AL: DFM100: 866199 Except product M1943AL: HeartStart MRx: 861288, 861289, 861464, 861482, 861491, 860396, 860397, 860398, 861465 | • | • |
| 989803105691 | M1943A | SpO ₂ Adapter Cable (8-pin), adapts 9-pin sensors to 8-pin sockets Length: 1 m (3.3 ft) | 1 cable | | • | • |
| 989803128651 | M1943AL | SpO ₂ Adapter Cable (8-pin), adapts 9-pin sensors to 8-pin sockets Length: 3 m (9.8 ft) | 1 cable | | • | • |
| 989803136591 | M1943NL | SpO ₂ Adapter Cable for OxiMax™ Instruments (8-pin), adapts 9-pin sensors to 8-pin sockets Use only with IntelliVue A02 and A04 sockets Length: 3 m (9.8 ft) | 1 cable | Pulse oximetry module: 862112; Nellcor OxiMax IV Module: 867193 MMS: 862442; MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MMX: 867036 | • | • |

Adapter cables

Adapter cables for Philips sensors



989803105681
M1941A

989803105691
M1943A

989803128651
M1943AL

989803136591
M1943NL

Adapter cables

Masimo® rainbow® SET™ SpO₂ adapter cables

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|--|--------|-----|
| 989803182411 | | rainbow® SET™ cable for rainbow, angled Length: 1.2 m (4 ft) | 1 cable | Masimo® rainbow® Module: 867191; MMS: 866462; X2: 865039; X3: 86703; MMX: 867036; MP2: 865040; MP5: 865024, 865322 | • | • |
| 989803186471 | | rainbow® SET™ cable for rainbow, straight Length: 3.7 m (12 ft) | 1 cable | Masimo® rainbow® Module: 867191; MMS: 866462; X2: 865039; X3: 86703; MMX: 867036; MP2: 865040; MP5: 865024, 865322; VS4: 863283 | • | • |
| 989803186501 | | rainbow® SET™ cable for LNCS® Length: 3 m (10 ft) | 1 cable | | • | • |
| 989803186531 | | rainbow® SET™ cable for LNCS®, straight Length: 1.2 m (4 ft) | 1 cable | MMS: 866462; MP2: 865040; MP5: 865024, 865322; X2: 865039; X3: 86703; MMX: 867036; Masimo® rainbow® Module: 867191; VS4: 863283; VS30: 863359, 863380; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304, 863322 | • | • |



989803182411



989803186471

Adapter cables

Masimo® rainbow® SET™ LNCS® SpO₂ adapter cables

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|---|--------|-----|
| 989803186501 | | rainbow® SET™ cable for LNCS® Length: 3 m (10 ft) | 1 cable | Masimo® rainbow® Module: 867191; MMS: 866462; X2: 865039; MMX: 867036; MP2: 865040; MP5: 865024, 865322; VS4: 863283 | • | • |
| 989803186531 | | rainbow® SET™ cable for LNCS®, straight Length: 1.2 m (4 ft) | 1 cable | MMS: 866462; MP2: 865040; MP5: 865024, 865322; X2: 865039; MMX: 867036; Masimo® rainbow® Module: 867191; VS4: 863283; VS30: 863359, 863380; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304, 863322 | • | • |



989803186501



989803186531

Adapter cables

Masimo® LNCS® SpO₂ adapter cable

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|--|--------|-----|
| 989803148221 | | Masimo® SET™ LNCS® adapter cable, adapts LNCS® sensors to Philips monitors Length: 3 m (10 ft) | 1 cable | VM1: 863264, 863265, 863266; VM4: 863063, 863085; VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088, 863317; VS2: 863278, 863279; VS3: 863069, 863070, 863071, 863072, 863073, 863074; VS4: 863283; VSi: 863276, 863277; MMS: 862442; MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MMX: 867036; MX40: 865350, 865351; Pulse oximetry module: 862112; Masimo® SET™ Module: 867192; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304, 863322; FM30: 862199; FM40: 865132; FM50: 865071 | • | • |



989803148221

Adapter cables

Masimo RD SET SpO₂ adapter cables

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-----------|---|--------|-----|
| 989803206901 | | RD SET MP-05 Philips Dual Keyed Cable; 5 ft. | 1 cable | MP2: 865040; MP5: 865024, 865322; MMS: 862442, 866462; MMX: 867036; MX100: 867033; X2: 865039; X3: 867030 | • | • |
| 989803206911 | | RD SET MP-12 Philips Dual Keyed Cable; 12 ft. | 1 cable | | • | • |

Adapter cables

Masimo RD rainbow SET SpO₂ adapter cables

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---------------------------------------|-----------|---|--------|-----|
| 989803206931 | | RD rainbow SET 25R – 12 Cable; 12 ft. | 1 cable | MP2: 865040; MP5: 865024, 865322; MMS: 862442, 866462; MMX: 867036; MX100: 867033; X2: 865039; X3: 867030 | • | • |



989803206901



989803206911



989803206931

Mobile CL non-invasive blood pressure cuffs

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|-----------------|--------|-----|
| 989803163201 | | Mobile CL single-patient use adult cuff, navy blue Cuff width: 13 cm (5.1 in), bladder length: 26 cm (10.2 in), limb circumference: 26–34.5 cm (10.2-13.5 in) | 20 cuffs | NBP pod: 865216 | • | • |



989803163201

Mobile CL non-invasive blood pressure accessories

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|------------|-----------------|--------|-----|
| 989803163131 | | Mobile CL extension air hose Length: 1 m (3.3 ft) | 1 hose | NBP pod: 865216 | • | • |
| 989803163251 | | Mobile CL NBP cradle kit | 20 cradles | | • | • |
| 989803187431 | | Mobile CL adult bayonet air hose Length: 1 m (3.3 ft) | 1 hose | | • | • |



989803163131



989803163251



989803187431

Temperature

Continuous measurement

Single use probes

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|--|--------|-----|
| 989803105321 | M1837A | Esophageal/rectal temperature probe, 9 FR, single-patient use, sterilized, continuous monitoring Length: 91 cm (36 in) | 20 probes | MP monitors: MP2: 865040; MP5: 865024; SureSigns monitors: VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; Temperature module: 862291; Microstream ext.: 862393, 865377 | | • |
| 989803100941 | 21090A | Esophageal/rectal temperature probe, 12 FR, single-patient use, sterilized, continuous monitoring Length: 91 cm (36 in) | 20 probes | | | • |



989803105321
M1837A



989803100941
21090A

Temperature

Continuous measurement

Reusable extension adapter cables for single-patient probes

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|----------|--------|-----|
| 989803100921 | 21082A | Adapter cable, connects single-use temperature probes to Philips (or HP or Agilent) monitors, accepts 3.5 mm (1/8 in) mini phone plug, 2-pin connector Length: 3 m (10 ft) | 1 cable | | | • |
| 989803100931 | 21082B | Adapter cable, connects single-use temperature probes to Philips (or HP or Agilent) monitors, accepts 3.5 mm (1/8 in) mini phone plug, 2-pin connector Length: 1.5 m (5 ft) | 1 cable | | | • |



989803100921
21082A



989803100931
21082B

Temperature

Continuous measurement

Reusable autoclavable probes

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|---|--------|-----|
| 989803203561 | | Esophageal/rectal temperature probe (autoclavable), 12 FR, reusable, continuous monitoring, 2-prong plug Length: 3 m (10 ft) | 1 probe | MP monitors: MP2: 865040; MP5: 865024; SureSigns monitors: VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088; CM10: 863301, 863310; CM12: 863303, 863311; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322; | • | • |
| 989803203571 | | Esophageal/rectal temperature probe (autoclavable), 10 FR, reusable, continuous monitoring, 2-prong plug Length: 3 m (10 ft) | 1 probe | MMS: 862442; X2: 865039; X3: 867030; MX100: 867033; Temperature module: 862291; Hemodynamic ext.: M3012A; 867039 Capnography ext.: M3014A; 867040 Microstream ext.: 862393, 865377, 867041 MMX: 867036 | • | • |
| 989803203581 | | Skin surface temperature probe (autoclavable), tape-on, reusable, continuous monitoring, 2-prong plug Length: 3 m (10 ft), diameter: 9.5 mm (0.37 in) | 1 probe | | • | • |



989803203561



989803203571



989803203581

Temperature

Spot check measurement

Predictive reusable temperature probes

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|-----------|---|--------|-----|
| 989803143381 | | Reusable oral temperature probe and well kit, predictive temp monitoring 2.7 m (9 ft) | 1 probe | VS2: 863278, 863081, 863082; VS3: 863070, 863073, 863074; VS4: 863276; VM4: 863063, VM6: 863064, 863065, VM8: 863066, MP5: 865024, 865322, 865120 | • | |
| 989803143391 | | Reusable rectal temperature probe and well kit, predictive temp monitoring 2.7 m (9 ft) | 1 probe | VS2: 863278, 863081, 863082; VS3: 863070, 863073, 863074; VM4: 863063, VM6: 863064, 863065, VM8: 863066, MP5: 865024, 865322, 865121 | • | |



989803143381



989803143391

Batteries

Patient monitoring batteries

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|-------------|--|--------|-----|
| 451261017451 | M4607-60001 | Battery charger adapter | 1 adapter | X2: 865039; MP2: 865040 | • | • |
| 989803135861 | M4605A | Rechargeable lithium-ion battery, 10.8 volts Shelf life: at least 6 months | 1 battery | IntelliVue MP20, MP30, MP40, MP50; MP5: 865024, 865322, 865120; FM20: 862198; FM30: 862199; Battery ext.: 865297; MX400: 866060; MX450: 866062; MX500: 866064; MX550: 866066 | • | • |
| 989803148701 | M4607A | Rechargeable lithium-ion battery, 10.8 volts Shelf life: at least 4 months | 1 battery | X2: 865039; MP2: 865040 | • | • |
| 989803174131 | | MX40 lithium-ion battery, 3-pack, 3.7 volts, rechargeable Shelf life: 3 months | 3 batteries | MX40: 865350, 865351, 867146 | • | • |
| 989803174891 | | MX40 battery adapter 3-pack, AA, disposable | 1 adapter | | • | • |
| 989803176201 | | MX40 Lithium ion battery 1-pack, 3.7 volts, rechargeable Shelf life: 3 months | 1 battery | | • | • |
| 989803174881 | | VM lithium-ion battery Shelf life: at least 6 months | 1 battery | VS: 863275, 863276, 863277; VS2: 863278, 863279; | • | • |
| 989803189981 | | Lithium-ion 3-cell battery pack Shelf life: at least 6 months | 1 battery | CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304 | • | |
| 989803194541 | | Lithium-ion 9-cell battery pack Shelf life: at least 6 months | 1 battery | CM10: 863301, 863365, 863310; CM12: 863303, 863366, 863311; CM100: 863300; CM120: 863302; CM150: 863304; Vital Sign View Station: 863067; VM4: 863063, 863085, 863287; VM6: 863064, 863065, 863086, 863288; VM8: 863066, 863068, 863087, 863088, 863317, 863289; VS2: 863278, 863279, 863079, 863080, 863081, 863082; VS3: 863069, 863070, 863071, 863072, 863073, 863074, 863285; VS4: 863283, 863286 | • | • |
| 989803196521 | | X3/MX100 lithium-ion battery Shelf life: at least 4 months | 1 battery | X3: 867030; MX100: 867033 | • | • |
| 989803202201 | | X3/MX100 battery charger adapter Shelf life: at least 4 months | 1 adapter | | • | • |
| 989803199221 | | Lithium Ion Battery Pack 9 Cells Shelflife: at least 6 months | 1 battery | 860354, 860355, 860429, 860437, 863359, 863380, 863300, 863301, 863302, 863303, 863304, 863320, 863365, 863321, 863366, 863322 | • | • |

Batteries

Patient monitoring batteries

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|---|--------------|---------------------------|--------|-----|
| 989803110231 | ULBU9VLJ | Lithium battery, 9 volts Shelf life: 10 months | 10 batteries | Telemetry: M2600A, M2601A | • | • |



451261017451
M4607-60001

989803135861
M4605A

989803148701
M4607A

989803174131

989803174881



989803174891

989803176201

989803194541

989803196521

989803110231
ULBU9VLJ

Accessories

Patient monitoring cases, bags and pads

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|-------------------------|-----------------|---|--------|-----|
| 989803174141 | | Carry pouch, waterproof | 50 pouches/box | MX40: 865350, 865351, 865352; NBP pod: 865216 | • | • |
| 989803174151 | | Carry pouch, waterproof | 200 pouches/box | | • | • |



989803174141 and 989803174151

Accessories

Telemetry pouches and protective covers

| Ordering code | Alternative ordering code | Description | Packaging | Use with | Europe | USA |
|---------------|---------------------------|--|--|---|--------|-----|
| 989803101971 | 9300-0768-050 | Disposable water and chemical resistant telemetry pouch with snaps | 50 pouches/box | Transmitter: 78100A, M1400A; Telemetry: M2600A; NBP pod: 865216 | • | • |
| 989803101981 | 9300-0768-200 | Disposable water and chemical resistant telemetry pouch with snaps | 50 pouches/box 4 boxes/case (200 pouches/case) | | • | • |
| 989803137831 | | Transceiver pouch shoulder strap and viewing window | 50 pouches/box | TRx4851A: 862231; NBP pod: 865216 | • | • |
| 989803140371 | | Transceiver pouch shoulder strap and viewing window | 50 pouches/box 4 boxes/case (200 pouches/case) | | • | • |



989803137831 and 989803140371



Patient monitors and defibrillators

Product list

Patient monitors

| Ordering code | Alternative ordering code | Description |
|--|---------------------------|---|
| 862108 | | IntelliVue M4841A telemetry transceiver |
| 862231 | | IntelliVue 2.4 GHz telemetry transceiver |
| 862318 | M1176A | Merlin bedside monitor |
| 862439 | | IntelliVue 1.4 GHz telemetry transceiver |
| 862442 | M3001A | IntelliVue multi-measurement server |
| 862474, 862478 | | C3 patient monitor |
| 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863062 | | SureSigns VS1 monitor |
| 863063, 863085, 863287 | | SureSigns VM4 bedside patient monitor |
| 863064, 863065, 863086, 863288 | | SureSigns VM6 bedside patient monitor |
| 863066, 863068, 863087, 863088, 863317, 863289 | | SureSigns VM8 bedside patient monitor |
| 863069, 863070, 863071, 863072, 863073, 863074, 863285 | | SureSigns VS3 vital signs monitor |
| 863077 | | SureSigns VM3 patient monitor |
| 863275, 863276, 863277 | | SureSigns VS1 vital signs monitor |
| 863278, 863279, 863079, 863080, 863081, 863082 | | SureSigns VS2 vital signs monitor |
| 863283, 863286 | | SureSigns VS4 vital signs monitor with NBP and SpO ₂ |
| 863300, 863320 | | Efficia CM100 patient monitor |
| 863301 | | Efficia CM10 patient monitor |
| 863365 | | Efficia CM10 patient monitor (Brazil) |
| 863302, 863321 | | Efficia CM120 patient monitor |
| 863303 | | Efficia CM12 patient monitor |
| 863366 | | Efficia CM12 patient monitor (Brazil) |
| 863304, 863322 | | Efficia CM150 patient monitor |

Patient monitors

| Ordering code | Alternative ordering code | Description |
|----------------|---------------------------|---|
| 863359, 863380 | | EarlyVue VS30 vital signs monitor |
| 865024 | M8105A | IntelliVue MP5 portable patient monitor |
| 865039 | M3002A | IntelliVue MMS X2 measurement module and patient monitor |
| 865040 | M8102A | IntelliVue MP2 portable patient monitor |
| 865120 | M8105AT | IntelliVue MP5T telemetry system and patient monitor |
| 865215 | | IntelliVue Cableless SpO ₂ pod |
| 865216 | | IntelliVue Cableless NBP pod |
| 865221 | | IntelliVue Cableless transmitter |
| 865322 | M8105AS | IntelliVue MP5SC spot check patient monitor |
| 865350 | | IntelliVue MX40 wearable patient monitor (US only), 1.4 GHz Smart-hopping |
| 865351, 867146 | | IntelliVue MX40 wearable patient monitor, 2.4 GHz Smart-hopping |
| 865352 | | IntelliVue MX40 wearable patient monitor, 802.11a/b/g |
| 866060 | | IntelliVue MX400 portable/bedside monitor |
| 866062 | | IntelliVue MX450 portable/bedside patient monitor |
| 866064 | | IntelliVue MX500 portable/bedside monitor |
| 866066 | | IntelliVue MX550 portable/bedside patient monitor |
| 866462 | M3001AL | IntelliVue multi-measurement server SLCP |
| 867030 | | IntelliVue X3 measurement module and patient monitor |
| 867033 | | IntelliVue MX100 patient monitor |
| 867041 | | Microstream Extension |

Measurement modules

| Ordering code | Alternative ordering code | Description |
|---------------|---------------------------|---|
| 453563372301 | M1025A | Anesthesia gas module |
| 453563372311 | M1025B | Anesthesia gas module |
| 989803110391 | M1001A | 3 ECG channel parameter module |
| 989803110421 | M1008A | Noninvasive blood pressure module |
| M1002A | M1002A | 3 ECG lead/respiration module |
| 860335 | | XPER Flex Cardio medical equipment |
| 860338 | | XPER Flex Cardio medical equipment |
| 862111 | M3012A | Hemodynamic measurement server extension |
| 862112 | M1020B | IntelliVue pulse oximetry module |
| 862119 | M1026B | Anesthesia gas module |
| 862165 | M1013A | IntelliVue G1 anesthesia gas module |
| 862187 | M3014A | Capnography extension |
| 862227 | M1019A | IntelliVue G5 anesthesia gas module |
| 862275 | M1001B | 3 ECG channel parameter module |
| 862276 | M1002B | 3 ECG channel/respiration module |
| 862278 | M1008B | Universal NBP module |
| 862279 | M1012A | Cardiac output module |
| 862281 | M1015A | Sidestream module |
| 862282 | M1016A | Mainstream CO ₂ module |
| 862283 | M1018A | Transcutaneous pO ₂ /pCO ₂ parameter module |
| 862286 | M1026A | Anesthesia gas module |
| 866078 | M1027B | IntelliVue EEG module |
| 862291 | M1029A | Temperature module |
| 862393 | M3015A | Microstream [®] CO ₂ extension |
| 862394 | M3016A | Mainstream CO ₂ 2nd press temp measurement extension |

Measurement modules

| Ordering code | Alternative ordering code | Description |
|---------------|---------------------------|--|
| 865021 | M1014A | IntelliVue spirometry module |
| 865377 | M3015B | Microstream [®] CO ₂ extension |
| 865383 | | IntelliVue NMT module |
| 866044 | | MB1 invasive and temp bundle |
| 866045 | | MB2 etCO ₂ bundle |
| 866046 | | MB3 basic bundle |
| 866149 | | Tympanic temperature module |
| 866173 | | IntelliVue G7m anesthesia gas module |
| 867039 | | IntelliVue hemodynamic extension |
| 867040 | | Capnography extension |
| 867191 | | Masimo rainbow [®] SET™ IV module |
| | M1006A, M1006B | Invasive blood pressure module |

Fetal monitors

| Ordering code | Alternative ordering code | Description |
|---------------|---------------------------|---|
| 862198 | M2702A | Avalon FM20 fetal monitor |
| 862199 | M2703A | Avalon FM30 fetal monitor |
| 862334 | M1350B | Series 50 XM fetal/maternal monitor |
| 862335 | M1350C | Series 50 XMO fetal/maternal monitor |
| 862336 | M1351A | Series 50A antepartum fetal monitor |
| 862337 | M1353A | Series 50IP-2 intrapartum fetal monitor |
| 865071 | M2705A | Avalon FM50 fetal monitor |
| 865132 | M2704A | Avalon FM40 fetal monitor |
| 867245 | | Avalon Toco+ MP transducer |
| 867246 | | Avalon US transducer |
| 867247 | | Avalon ECG/IUP transducer |
| 867248 | | Avalon TOCO MP transducer |
| 867249 | | Avalon Toco+ transducer |
| 866075 | | Avalon Cableless Toco+ MP transducer |
| 866076 | | Avalon Cableless ultrasound transducer |
| 866077 | | Avalon Cableless ECG/IUP transducer |

Defibrillators

| Ordering code | Alternative ordering code | Description |
|--------------------------------|---------------------------|--|
| 860396 | M3536M7 | HeartStart MRx pre-hospital government without temperature |
| 860397 | M3536M8 | HeartStart MRx pre-hospital government with temperature |
| 860398 | M3536M9 | HeartStart MRx pre-hospital government critical care |
| 861266 | M4735A | HeartStart XL defibrillator/monitor |
| 861288, 861493, 861494, 861495 | M3535A | HeartStart MRx hospital monitor/defibrillator |
| 861289 | M3536A | HeartStart MRx pre-hospital monitor/defibrillator |
| 861290 | | HeartStart XL+ defibrillator/monitor |
| 861464 | M3536M | HeartStart MRx monitor/defibrillator, airworthy no case |
| 861465 | M3536MC | HeartStart MRx monitor/defibrillator, airworthy with case |
| 861480 | M3536M1 | HeartStart MRx crash cart monitor/defibrillator |
| 861481 | M3536M2 | HeartStart MRx ICU monitor/defibrillator |
| 861482 | M3536M3 | HeartStart MRx EMS responder monitor/defibrillator |
| 861483 | M3536M4 | HeartStart MRx EMS responder monitor/defibrillator |
| 861484 | M3536M5 | HeartStart MRx EMS CC monitor/defibrillator |
| 861491 | M3536M6 | HeartStart MRx EMS CCT monitor/defibrillator, no Bluetooth |
| 866199 | | Efficia DFM100 defibrillator/monitor |
| 867172 | | Heartstart Intrepid monitor/defibrillator |

